



# Gentle to animals

Good for animals,  
made for the practice.



# Treating animals well



As an innovator with decades of experience in human pre-analytics, SARSTEDT offers a number of solutions that have also been proven in veterinary medicine and where the focus is on hygiene, user friendliness and animal welfare.

With products such as the S-Monovette®, we set hygiene standards for closed venous blood collection from pets such as dogs and cats. We also offer suitable solutions for collecting samples that meet the specific requirements of the smallest patients such as small mammals, birds or reptiles.

Test our well conceived systems and see for yourself the many practical benefits for your day-to-day work.

Blood collection with SARSTEDT. Good for animals, made for the practice.




“We use the S-Monovette® for large dogs. This allows us to fill the sample tube easily and cleanly even with just one hand.”

## Pre-analytics made by SARSTEDT

Pre-analytics is the most decisive phase for the examination of every animal, as even with the highest precision of modern analysis equipment, the result is only as good as the quality of the sample.

With our expertise in pre-analytics, we assist you and give you products that support your work routine. The goal is to give you more time for the animals and owners that come to your practice.



“If I get only a few drops of blood, the Micro needle helps me get as much as possible from the collection – it prevents coagulation at the needle aperture.”



“Using the Multifly® needle is an enormous relief for me. When the puncture is made, I have enough time to fill the required sample tubes.”



“We call the Multivette® the ‘vacuum’. We just connect it to the Luer needle and can fill two or three tubes with blood while we are in the vein.”

- ✓ Safe
- ✓ Gentle
- ✓ Hygienic

The convenient solution for every animal – together for greater animal welfare.

Many different patients come to your veterinary practice, all with their individual needs. SARSTEDT gives you suitable solutions for pre-analytics in veterinary medicine.

- The S-Monovette® – safe, gentle and hygienic blood collection from larger pets
- The Multivette® – collect blood gently and easily from small pets using only venous pressure
- The Micro needle – even in especially sensitive animals, collect every drop reliably, e.g. with a Micro sample tube or a Microvette®
- Many other practical systems and single products for collecting blood, urine and faeces samples

More about Animal Health Care from SARSTEDT



animalhealthcare.sarstedt.com



# Closed blood collection from larger pets




The S-Monovette® blood collection system allows you to perform **closed blood collection with two techniques**. The aspiration technique enables gentle collection adapted to the blood flow. When blood flow is strong, the dual system of the S-Monovette® also allows the conventional vacuum technique for blood collection.

## The sustainable blood collection system

The great advantage of blood collection with an S-Monovette® is that the blood can be taken directly into the correctly prepared sample tube. No extra disposable syringe is needed to transfer the blood sample to other tubes. This saves resources and makes blood collection safer and more hygienic than alternative collection options.

The more sustainable cardboard packaging also keeps disposal costs and waste volume low.

Make the trip to the veterinarian a positive experience for your patients and their owners



Test now  
for free!



animalhealthcare.  
sarstedt.com




*"I went to the doc again today, but this time something was different: Everything was completely relaxed.*

*Before, I just got pricked in the leg and this red liquid spilled out everywhere. Today, everything went very quickly without any mess. That must be because of the doc's new equipment.*

*I hope it's like this every time now!"*

**Lucky, male boxer**



## The S-Monovette®

Safety needle guard for use with one hand

Secure, tight connection due to the pre-mounted adapter

Visible check for successful vein puncture using the first drop of blood that enters the S-Monovette® or the tubing (with the Safety-Multifly® needle).

The S-Monovette® is available for blood volumes from 1.2–4.5 ml in different preparations.

Gentle collection using the aspiration technique.



**By the way** the S-Monovette® can also be used with the Safety-Multifly® needle.

# Relaxed blood collection from small pets




The Multivette® 600 is designed and prepared **for small blood volumes** of 600 µl. Blood is collected in a nearly closed condition and utilises natural venous pressure. This makes collection especially gentle and easy.

The blood sample can be centrifuged directly in the Multivette® 600. This makes blood collection easier and saves time.

The small inner diameter makes it particularly easy to pipette blood after centrifugation. Additionally, the sample can be securely closed for sending to the laboratory.


## The Multivette® 600



*"Today my can openers\* put me in this box again and took me away. When they finally let me out, they first rubbed me behind my ears.*

*Then something magical happened: They held this thing against my fur, the colour changed from white to red, they rubbed me again and I got to go home. It wasn't bad at all."*

**Bella, Bengal cat**  
*\*In the cat world, can opener means "human"*




The inner tube with a special design is filled from venous pressure alone.

Different preparations for 600 µl blood samples



# Successful drip method in especially sensitive animals



"Hi guys! My mom was worried because I haven't been feeling well for a few days. My vet first drew some blood to find out why. She had a special needle for me that I could hardly feel. Actually, I didn't really notice anything.

Mom was also much more relaxed than usual. That was most certainly because it was over so quickly – she can't stand the sight of blood..."

**Coco, female Chihuahua**



## The Micro needle with Micro sample tubes or Microvette®

Especially in very small animals and **difficult vein conditions**, every drop of blood counts. The special Micro needle ensures that every drop of blood flows into the sample tube and does not coagulate beforehand. The appropriate sample tube can be selected according to the expected sample volume.

The Micro sample tubes are suitable for volumes of 1.3 ml. For smaller patients, the Microvette® in volumes 100 – 500 µl is the right choice.

Optimal flow properties

Minimal blood loss with no residual volume in the attachment

Direct filling avoids blood loss through later transfer

Flexibility thanks to different volumes and preparations



Test now for free!



animalhealthcare.  
sarstedt.com



**SARSTEDT**

## SARSTEDT AG & Co. KG

P.O. Box 12 20  
D-51582 Nümbrecht

Phone: +49 2293 305 0  
Fax: +49 2293 305 3992

Customer Service Germany  
Phone +49 0800 0 83 305 0

[export@sarstedt.com](mailto:export@sarstedt.com)  
[www.sarstedt.com](http://www.sarstedt.com)



Further information and free samples  
for testing:  
[animalhealthcare.sarstedt.com](http://animalhealthcare.sarstedt.com)

If you have any questions:  
We'd be happy to help!

Visit our website: [www.sarstedt.com](http://www.sarstedt.com)