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Intended use

Microvette® Capillary Blood Collection Systems are blood collection systems and are used with lancets to collect capillary blood using End-to-End capillaries or the collection rim. The Microvette® Capillary Blood Collection Systems are used for collecting, manual in-house transport and processing of capillary blood specimens in the clinical laboratory. The products are intended for use in a professional environment by qualified medical and laboratory personnel.

Product description

The Microvette® Capillary Blood Collection Systems consist of a plastic receptacle and a colour-coded plastic screw cap and various additives (preparations). The selection of the additive depends on the analytical test method. It is prespecified by the manufacturer of the test reagents and/or the analysis equipment on which the test is conducted.

Microvette® screw cap colour codes:

Additive	Letter code	Cap colour based on ISO 6710	Cap colour based on BS 4851*
Tube without preparation			
Microvette® 500 neutral Z in a cylindrical inner receptacle	Z	red	-
Serum tube			
Microvette® 100 with serum in a conical inner receptacle	CAT	red	white
Microvette® 200 with serum in a conical inner receptacle	CAT	red	white
Microvette® 200 with serum gel in a cylindrical inner receptacle	CAT	brown	brown
Microvette® 300 with serum in a conical inner receptacle	CAT	red	white
Microvette® 500 with serum in a cylindrical inner receptacle	CAT	red	white
Microvette® 500 with serum gel in a cylindrical inner receptacle	CAT	brown	brown
Microvette® 500 with serum gel in a cylindrical inner receptacle**	CAT	white	-
Lithium heparin tube			
Microvette® 100 with lithium heparin in a conical inner receptacle	LH	green	orange
Microvette® 200 with lithium heparin in a conical inner receptacle	LH	green	orange
Microvette® 300 with lithium heparin in a conical inner receptacle	LH	green	orange
Microvette® 500 with lithium heparin in a cylindrical inner receptacle	LH	green	orange
Microvette® 500 with lithium heparin gel in a cylindrical inner receptacle	LH	green	orange
EDTA tube			
Microvette® 100 with K3 EDTA in a cylindrical inner receptacle	K3E	violet	red
Microvette® 200 with K3 EDTA in a cylindrical inner receptacle	K3E	violet	red
Microvette® 200 with K2 EDTA in a cylindrical inner receptacle	K2E	violet	-
Microvette® 500 with K3 EDTA in a cylindrical inner receptacle	K3E	violet	red
Microvette® 500 with K2 EDTA in a cylindrical inner receptacle	K2E	violet	-
Tube with glycolysis inhibitors			
Microvette® 200 with fluoride/heparin in a conical inner receptacle	FH	grey	yellow
Microvette® 300 with fluoride/heparin in a conical inner receptacle	FH	grey	yellow
Microvette® 500 with fluoride/heparin in a conical inner receptacle	FH	grey	-

*British Standard BS4851:1982 Specification for single use labelled medical specimen containers for haematology and biochemistry.
 **brown outer container with light protection properties

SARSTEDT Microvette® neutral

The Microvette® neutral can be used for user-specific tests.

SARSTEDT Microvette® Serum/Serum Gel

The Microvette® with serum contains a clotting activator and is used for collecting serum. The Microvette® with serum gel contains a clotting activator and is used for collecting serum. After centrifugation, the gel forms a layer separating the coagulum and the serum. Microvettes with serum / serum gel are intended for clinical chemistry tests from serum.

SARSTEDT Microvette® lithium heparin / lithium heparin gel

The Microvette® with lithium heparin is coated with lithium heparin as an anticoagulant and is used for collecting plasma. The Microvette® with lithium heparin gel contains lithium heparin as an anticoagulant and is used for collecting plasma. After centrifugation, the gel forms a layer separating the plasma and the cellular blood components. Microvettes with lithium heparin / lithium heparin gel are intended for clinical chemistry tests from plasma.

SARSTEDT Microvette® K3/K2 EDTA

The Microvette® with K3/K2 EDTA is coated with the anticoagulant K3 EDTA/K2 EDTA and is intended for haematological tests.

SARSTEDT Microvette® Fluoride/heparin

The Microvette® with fluoride/heparin is coated with the anticoagulant heparin and contains fluoride for glycolysis inhibition. It is used to test glucose levels.

NOTE: Fluoride (used for measuring glucose) can cause an increase in haemolysis. Further information on substances that can cause interference can be found in the respective instructions for use of the assay manufacturer.

Safety information and warnings

1. General precautions: Use gloves and other general personal protective equipment to protect yourself against blood and possible exposure to pathogens transmitted by blood.
2. Handle all biological specimens and sharp blood collection utensils according to the guidelines and procedures in your facility. In case of direct contact with biological specimens or a needlestick injury, consult a doctor due to the risk of transmission of HIV, HCV, or HBV or other infectious diseases. Always follow the safety guidelines and procedures of your facility.
3. Dispose of all sharp objects for capillary blood collection in appropriate disposal containers.
4. Under- or overfilling the Microvettes leads to an incorrect ratio of blood to the preparation/additive and can result in incorrect analysis results.
5. The blood collected and processed with the Microvette® is not intended to be re-injected into the human body.
6. Do not use the product after the expiry date. The expiry date of the product is the last day of the month and year as indicated.

Storage

Store the products at room temperature.

Limitations

1. If blood samples are stored in the Microvette®, the stability of the analytes should be evaluated by the respective laboratory or taken from the instructions for use of the analysis device manufacturer.
2. Even if plasma or serum is separated by centrifuging the Microvettes and/or a barrier is present, all cells are not necessarily completely separated. Residual metabolism or natural degradation can affect the concentrations of analytes. The stability of the analyte should be evaluated with respect to the storage container and the conditions of the respective laboratory.
3. For therapeutic medication, the suitability of the sample material must be checked in the manufacturer's instructions for use of the assay / analysis device.

Collecting and handling specimens

READ THIS DOCUMENT COMPLETELY BEFORE STARTING TO DRAW CAPILLARY BLOOD.

Preparation for drawing capillary blood and material needed:

1. All necessary Microvettes, labelled by size and additive.
2. Gloves, gown, goggles or other suitable protective wear to protect against bloodborne pathogens or potentially infectious materials.
3. Labels for specimen identification.
4. (Safety) lancet or incision lancet.
5. Disinfectant for cleaning the puncture site (follow the facility's guidelines for preparing the puncture site for taking blood samples). Do not use any alcohol-based cleaning materials if the specimens are to be used for a blood alcohol test.
6. Dry, clean swabs.
7. Plasters.
8. Sharps disposal container for the safe disposal of used material.

Recommended order of collection (CLSI* GP42):

1. Tubes with EDTA
2. Tubes with lithium heparin with/without gel
3. Tubes with glycolysis inhibitors
4. Tubes for coagulation tests of serum with/without gel

NOTE: The first drop of blood is discarded.


NOTE: Follow the regulations of your facility for the order of collection.

Capillary technique with the End-to-End capillary

General instructions:

WEAR GLOVES TO MINIMISE THE RISK OF EXPOSURE DURING CAPILLARY BLOOD COLLECTION.

1. Select a suitable Microvette® for the required specimen.
2. Select a site for the puncture with a (safety) lancet or incision lancet. Please follow the guidelines of the facility for preparing the puncture site.

- 1 Hold the Microvette® horizontally or tilted slightly and collect the drops of blood with the 100 µl or 200 µl End-to-End capillary.
- 2 Stop collecting blood when the End-to-End capillary is completely filled with blood.
- 3 Hold the Microvette® upright so the blood flows into the collection container.
- 4 Twist slightly to remove the cap incl. the End-to-End capillary and discard as one unit.
- 5 Remove the cap attached to the bottom of the tube and seal the tube ("click" position).
- 6 Invert the specimen carefully.

Information on coagulation / on the Microvette® Serum:


After collecting blood, let the the Microvette® Serum/Serum Gel coagulate for 30 minutes and the Microvette® neutral for at least 60 minutes in an upright position to ensure a clear separation layer after centrifugation. The recommended times are based on an intact coagulation process. Blood from patients with coagulation disorders due to their illness or under anticoagulant therapy requires more time for complete coagulation.

Blood collection with the collection rim

General instructions:

WEAR GLOVES TO MINIMISE THE RISK OF EXPOSURE DURING BLOOD COLLECTION.

1. Select a suitable Microvette® for the required specimen.
2. Select a site for the puncture with a (safety) lancet or incision lancet. Please follow the guidelines of the facility for preparing the puncture site.

- 1 Twist the cap slightly to take it off.
- 2 Attach the cap to the bottom of the tube.
- 3 Collect the drops of blood using the collection rim.
- 4 Remove the cap attached to the bottom of the tube and seal the Microvette® ("click" position).
- 5 Invert the specimen carefully.

Information on coagulation / on the Microvette® Serum:

After collecting blood, let the Microvette® Serum / Serum Gel coagulate for 30 minutes and the Microvette® neutral for at least 60 minutes in an upright position to ensure a clear separation layer after centrifugation. The recommended times are based on an intact coagulation process. Blood from patients with coagulation disorders due to their illness or under anticoagulant therapy requires more time for complete coagulation.

Centrifugation

WARNING!
Centrifuging Microvettes with cracks or centrifuging at an excessive centrifugal force can lead to breakage of the Microvettes, which can release potentially hazardous substances.

Please see our centrifugation recommendations for the various preparations in the table below.
Centrifuge inserts are to be selected according to the size of the Microvettes used. The relative centrifugal force has the following ratio to the selected rpm:

$$RCF (G-Force) = 11.2 \times r \times (RPM/1000)^2,$$

"RCF": "Relative centrifugal force"

"RPM": "Revolutions per minute"

"r" [in cm]: "Radius of rotation from the middle of the centrifuge to the bottom of the Microvette®,"

Microvettes without gel can be centrifuged in centrifuges with a fixed-angle or swing-out rotor.

Microvettes with gel are designed solely for use in centrifuges with swing-out rotor. Centrifugation in centrifuges with a fixed-angle rotor is not validated by SARSTEDT and is not recommended. In addition, centrifugation should be conducted no later than 2 hours after collection, as this could otherwise yield false results.

Microvettes should be centrifuged in accordance with the centrifugation conditions listed below. If other conditions are used, they must be validated by the user.

It must be ensured that the Microvettes fit well in the centrifuge inserts. Microvettes that protrude beyond the insert can catch on the centrifuge head and be broken. Ensure that the centrifuge is filled evenly. Please refer to the centrifuge instructions for use for further information.

CAUTION! Do not remove broken Microvettes by hand.

Refer to the centrifuge instructions for use for information on disinfecting the centrifuge.

Image	Preparation	Times and g values (=RCF*)				
		Standard recommendation		Alternative recommendation		
	Serum	5 min	10,000 x g	or	10 min	
	Lithium heparin					
	EDTA ²		2,000 x g			
	Glucose fluoride					
	Serum gel ¹		10,000 x g			4,000 - 10,000 x g
	Lithium heparin gel ¹					

¹ For tubes prepared with gel, we recommend using swing-out rotors.

² Only for plasma separation, not for haematological tests.

*RCF: Relative centrifugal force

Centrifugation at 20 °C

Disposal

1. The general hygiene guidelines and regulations for the proper disposal of infectious material must be observed and complied with.
2. Disposable gloves prevent the risk of infection.
3. Contaminated or filled capillary blood collection systems must be disposed of in suitable containers for hazardous biological waste that can be subsequently autoclaved and incinerated.
4. They must be disposed of in a suitable incinerator or by autoclaving (steam sterilisation).

Product-specific standards and guidelines as amended/updated

CLSI* GP39 "Tubes and Additives for Venous and Capillary Blood Specimen Collection", Approved Standard.

CLSI* GP42 "Collection of Capillary Blood Specimens", Approved Standard.

CLSI* GP44 "Procedures for the Handling and Processing of Blood Specimens for Common Laboratory Tests", Approved Guideline.

WHO/DIL/LAB/99.1 Rev.02 "WORLD HEALTH ORGANIZATION, et al. Use of anticoagulants in diagnostic laboratory investigations. Geneva: World Health Organization, 2002".

*CLSI (Clinical and Laboratory Standards Institute)

Key for symbols and labels:

	Article number
	Batch number
	Use by
	CE marking
	<i>In vitro</i> diagnostic device
	Follow the instructions for use
	If reused: Risk of contamination
	Keep away from sunlight
	Store in a dry place
	Manufacturer
	Country of manufacture

Technical modifications reserved.

All serious incidents relating to the product must be reported to the manufacturer and the competent national authority.

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CLSI* GP39 "Tubes and Additives for Venous and Capillary Blood Specimen Collection", Approved Standard.
 CLSI* GP42 "Collection of Capillary Blood Specimens", Approved Standard.
 CLSI* GP44 "Procedures for the Handling and Processing of Blood Specimens for Common Laboratory Tests", Approved Guideline.
 WHO/DIL/LAB/99.1 Rev02 "WORLD HEALTH ORGANIZATION, et al. Use of anticoagulants in diagnostic laboratory investigations. Geneva: World Health Organization, 2002".
 *CLSI (Clinical and Laboratory Standards Institute)

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- CLSI* GP39 "Tubes and Additives for Venous and Capillary Blood Specimen Collection", Approved Standard.
- CLSI* GP42 "Collection of Capillary Blood Specimens", Approved Standard.
- CLSI* GP44 "Procedures for the Handling and Processing of Blood Specimens for Common Laboratory Tests", Approved Guideline.
- WHO/DIL/LAB/99.1 Rev02 "WORLD HEALTH ORGANIZATION, et al. Use of anticoagulants in diagnostic laboratory investigations. Geneva: World Health Organization, 2002".
- *CLSI (Clinical and Laboratory Standards Institute)

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