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## SARSTEDT Microvette® APT K2E

<b>DE</b>	Gebrauchshinweis – SARSTEDT Microvette® APT K2E	2 - 6
<b>EN</b>	Instructions for use – SARSTEDT Microvette® APT K2E	7 - 11
<b>BG</b>	Инструкции за употреба – SARSTEDT Microvette® APT K2E	12 - 16
<b>CS</b>	Návod k obsluze – SARSTEDT Microvette® APT K2E	17 - 21
<b>DA</b>	Brugsanvisning – SARSTEDT Microvette® APT K2E	22 - 26
<b>EL</b>	Οδηγίες χρήσεως – SARSTEDT Microvette® APT K2E	27 - 31
<b>ES</b>	Instrucciones de uso – SARSTEDT Microvette® APT K2E	32 - 36
<b>ET</b>	Kasutusjuhend – Verevõtusüsteem SARSTEDT Microvette® APT K2E	37 - 41
<b>FR</b>	Mode d'emploi – Microvette® APT K2E SARSTEDT	42 - 46
<b>HR</b>	Uputa za upotrebu – SARSTEDT Microvette® APT K2E	47 - 51
<b>HU</b>	Használati utasítás – SARSTEDT Microvette® APT K2E	52 - 56
<b>IT</b>	Istruzioni d'uso – Microvette® APT K2E SARSTEDT	57 - 61
<b>KO</b>	사용 설명서 – SARSTEDT Microvette® APT K2E	62 - 66
<b>LT</b>	Naudojimo instrukcijos – SARSTEDT „Microvette® APT K2E“	67 - 71
<b>LV</b>	Lietošanas norādes – SARSTEDT Microvette® APT K2E	72 - 76
<b>NL</b>	Gebruiksaanwijzing – SARSTEDT Microvette® APT K2E	77 - 81
<b>NO</b>	Bruksanvisning – SARSTEDT Microvette® APT K2E	82 - 86
<b>PL</b>	Instrukcja obsługi – Probówka Microvette® APT K2E SARSTEDT	87 - 91
<b>PT</b>	Instruções de utilização – SARSTEDT Microvette® APT K2E	92 - 96
<b>RO</b>	Instrucțiuni de utilizare – Microvette® APT K2E SARSTEDT	97 - 101
<b>RU</b>	Инструкция по применению – SARSTEDT Microvette® APT K2E	102 - 106
<b>SK</b>	Návod na Použitie – SARSTEDT Microvette® APT K2E	107 - 111
<b>SL</b>	Navodila za uporabo – Epruveta Microvette® APT K2E SARSTEDT	112 - 116
<b>SV</b>	Bruksanvisning – SARSTEDT Microvette® APT K2E	117 - 121
<b>TH</b>	คำแนะนำในการใช้งาน – ระบบการเจาะเลือดจากเส้นเลือดฝอย Microvette® APT K2E ของ SARSTEDT	122 - 126
<b>TR</b>	Kullanım için talimatlar – SARSTEDT Microvette® APT K2E	127 - 131
<b>ZH</b>	使用说明 – SARSTEDT Microvette® APT K2E	132 - 136

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## Intended use

The Microvette® APT K2E is a blood collection system and is used with lancets to collect capillary blood using an End-to-End capillary or a collection rim. The Microvette® APT K2E is used for the collection, anticoagulation with K2-EDTA, processing and transport of capillary blood samples in the clinical laboratory. The blood samples obtained with the Microvette® APT K2E can be used for both manual and automated haematological tests. The product is intended for use in a professional environment by qualified medical and laboratory personnel.

## Product description

The Microvette® APT K2E consists of a plastic receptacle and a colour-coded plastic cap with a pierceable membrane. It is available in two sizes, 250 and 500 µl. The Microvette® APT 250 K2E has an End-to-End capillary.

## Microvette® screw cap colour codes:

Additive	Letter code	Cap colour based on ISO 6710	Cap colour based on BS 4851*
<b>EDTA tube</b>			
Microvette® APT 250 K2E	K2E	Violet	Red
Microvette® APT 500 K2E	K2E	Violet	Red

\*British Standard BS4851:1982 Specification for single use labelled medical specimen containers for haematology and biochemistry.

## Safety information and warnings

1. General precautions: Use gloves and other general personal protective equipment to protect yourself against blood and possible exposure to pathogens transmitted by blood.
2. Handle all biological specimens and sharp blood collection utensils according to the guidelines and procedures in your facility. In case of direct contact with biological specimens or a needlestick injury, consult a doctor due to the risk of transmission of HIV, HCV, or HBV or other infectious diseases. Always follow the safety guidelines and procedures of your facility.
3. Dispose of all sharp objects for capillary blood collection in appropriate containers for sharps disposal.
4. Under- or overfilling the Microvettes leads to an incorrect ratio of blood to the preparation/additive and can result in incorrect analysis results.
5. The blood collected and processed with the Microvette® is not intended to be re-injected into the human body.
6. Do not use a Microvette® after the expiry date. The expiry date of a Microvette® is the last day of the month and year indicated.

## Storage

Store the product at room temperature.

## Transport

The product is used as a primary container under ADR P650.

## Limitations

1. If blood samples are stored in the Microvette® APT K2E, the stability of the analytes should be evaluated by the respective laboratory or taken from the instructions for use of the analysis device manufacturer.
2. For therapeutic medication, the suitability of the sample material must be checked in the manufacturer's instructions for use of the assay/analysis device.
3. When using the Microvette® APT K2E on automated analysis devices, always observe the information provided by the device manufacturer.

## Collecting and handling specimens

**READ THIS DOCUMENT COMPLETELY BEFORE STARTING TO DRAW CAPILLARY BLOOD.**

### Preparation for drawing capillary blood and material needed:

1. All necessary Microvettes, labelled by size and additive.
2. Gloves, gown, goggles or other suitable protective wear to protect against bloodborne pathogens or potentially infectious materials.
3. Labels for specimen identification.
4. (Safety) lancet.
5. Disinfectant to clean the collection site (follow guidelines of the facility for sterile specimen collection for preparing the collection site). Do not use any alcohol-based cleaning materials if the specimens are to be used for a blood alcohol test.
6. Dry, clean swabs.
7. Plasters.
8. Sharps disposal container for the safe disposal of used material.

## Recommended order of collection (CLSI\* GP42):

1. Tubes with EDTA
2. Tubes with lithium heparin with/without gel
3. Tubes with glycolysis inhibitors
4. Tubes for coagulation tests of serum with/without gel

**NOTE: The first drop of blood is discarded.**


**NOTE: Follow the regulations of your facility for the order of collection.**

## Capillary technique with the End-to-End capillary

### General instructions:

WEAR GLOVES TO MINIMISE THE RISK OF EXPOSURE DURING CAPILLARY BLOOD COLLECTION.

1. Select a suitable Microvette® APT 250 K2E for the required specimen.
2. Select the location for the puncture with a (safety) lancet. Please follow the guidelines of the facility for preparing the puncture site.

### A

- 1 Hold the Microvette® APT 250 K2E horizontally or tilted slightly and collect the drops of blood with the 250 µl End-to-End capillary.
- 2 Stop collecting blood when the End-to-End capillary is completely filled with blood.
- 3 Hold the Microvette® APT 250 K2E upright so the blood flows into the collection container.
- 4 Twist slightly to remove the cap incl. the End-to-End capillary and discard as one unit.
- 5 Close with the included cap.
- 6 The Microvette® APT 250 K2E is properly closed when you hear a click.
- 7 Gently shake the Microvette® APT 250 K2E from side to side.

### B

- 1 For the drop technique using the collection rim, twist slightly to remove the cap incl. the End-to-End capillary and discard as one unit.
- 2 Collect the drops of blood using the collection rim.
- 3 Close with the included cap.
- 4 The Microvette® APT 250 K2E is properly closed when you hear a click.
- 5 Gently shake the Microvette® APT 250 K2E from side to side.


## Blood collection with the drip rim

### General instructions:

WEAR GLOVES TO MINIMISE THE RISK OF EXPOSURE DURING BLOOD COLLECTION.

1. Select a suitable Microvette® APT 500 K2E for the required specimen.
2. Select the location for the puncture with a (safety) lancet. Please follow the guidelines of the facility for preparing the puncture site.

- 1 Remove the pre-installed cap of the Microvette® APT 500 K2E with a slight twisting movement.
- 2 Collect the drops of blood using the collection rim.
- 3 Reclose it with the cap.
- 4 The Microvette® APT 500 K2E is properly closed when you hear a click.
- 5 Gently shake the Microvette® APT 500 K2E from side to side.

## Centrifugation

**WARNING!**  
Centrifuging Microvettes with cracks or centrifuging at an excessive centrifugal force can lead to breakage of the Microvettes, which can release potentially hazardous substances.

Centrifuge inserts are to be selected according to the size of the Microvettes used. The relative centrifugal force has the following relationship to the set RPM (revolutions per minute):

$$RCF \text{ (g force)} = 11.2 \times r \times (\text{RPM}/1000)^2,$$

"RCF": "relative centrifugal force",

"RPM": "revolutions per minute" ( $\text{min}^{-1}$ ),

"r" [in cm]: "Radius of rotation from the middle of the centrifuge to the bottom of the Microvette®",

Microvettes without gel can be centrifuged in centrifuges with a fixed-angle or swing-out rotor.

Centrifugation in centrifuges with a fixed-angle rotor is not validated by SARSTEDT and is not recommended. In addition, centrifugation should be conducted no later than 2 hours after collection, as this could otherwise yield false results.

Microvettes should be centrifuged in accordance with the centrifugation conditions listed below. If other conditions are used, they must be validated by the user.

It must be ensured that the Microvettes fit well in the centrifuge inserts. Microvettes that protrude beyond the insert can catch on the centrifuge head and be broken. Ensure that the centrifuge is filled evenly. Please refer to the centrifuge instructions for use for further information.

**CAUTION!** Do not remove broken Microvettes by hand.

Refer to the centrifuge instructions for use for information on disinfecting the centrifuge.

Image	Preparation	Times and g values (=RCF*)				
		Standard recommendation		Alternative recommendation		
	K2 EDTA	9 min	2,000 x g	or	5 min	4,000 x g

\*RCF: Relative centrifugal force

Centrifugation at 20°C

## Disposal

1. The general hygiene guidelines and regulations for the proper disposal of infectious material must be observed and complied with.
2. Disposable gloves prevent the risk of infection.
3. Contaminated or filled capillary blood collection systems must be disposed of in suitable containers for hazardous biological waste that can be subsequently autoclaved and incinerated.
4. They must be disposed of in a suitable incinerator or by autoclaving (steam sterilisation).

## Product-specific standards and guidelines as amended/updated

CLSI\* GP39 "Tubes and Additives for Venous and Capillary Blood Specimen Collection", Approved Standard.

CLSI\* GP42 "Collection of Capillary Blood Specimens", Approved Standard.

CLSI\* GP44 "Procedures for the Handling and Processing of Blood Specimens for Common Laboratory Tests", Approved Guideline.

WHO/DIL/LAB/99.1 Rev02 "WORLD HEALTH ORGANIZATION, et al. Use of anticoagulants in diagnostic laboratory investigations. Geneva: World Health Organization, 2002".

\*CLSI (Clinical and Laboratory Standards Institute)

## Key for symbols and labels:

	Article number
	Batch number
	Use by
	CE marking
	<i>In vitro</i> diagnostic device
	Follow the instructions for use
	If reused: Risk of contamination
	Store away from sunlight
	Store in a dry place
	Manufacturer
	Country of manufacture

Technical modifications reserved.

All serious incidents relating to the product must be reported to the manufacturer and the competent national authority.

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CLSI\* GP39 „Tubes and Additives for Venous and Capillary Blood Specimen Collection”, Approved Standard.

CLSI\* GP42 „Collection of Capillary Blood Specimens”, Approved Standard.

CLSI\* GP44 „Procedures for the Handling and Processing of Blood Specimens for Common Laboratory Tests”, Approved Guideline.

WHO/DIL/LAB/99.1 Rev02 „WORLD HEALTH ORGANIZATION, et al. Use of anticoagulants in diagnostic laboratory investigations. Geneva: World Health Organization, 2002”.

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