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## SARSTEDT Multivette® 600

<b>DE</b>	Gebrauchshinweis – SARSTEDT Multivette® 600	2 - 7
<b>EN</b>	Instructions for use – SARSTEDT Multivette® 600	8 - 13
<b>BG</b>	Инструкции за употреба – SARSTEDT Multivette® 600	14 - 19
<b>CS</b>	Návod k obsluze – SARSTEDT Multivette® 600	20 - 25
<b>DA</b>	Brugsanvisning – SARSTEDT Multivette® 600	26 - 31
<b>EL</b>	Οδηγίες χρήσεως – SARSTEDT Multivette® 600	32 - 37
<b>ES</b>	Instrucciones de uso – SARSTEDT Multivette® 600	38 - 43
<b>ET</b>	Kasutusjuhend – SARSTEDT Multivette® 600	44 - 49
<b>FR</b>	Mode d'emploi – Multivette® 600 SARSTEDT	50 - 55
<b>HR</b>	Uputa za upotrebu – SARSTEDT Multivette® 600	56 - 61
<b>HU</b>	Használati utasítás – SARSTEDT Multivette® 600	62 - 67
<b>IT</b>	Istruzioni d'uso – Multivette® 600 SARSTEDT	68 - 73
<b>KO</b>	사용 설명서 – SARSTEDT Multivette® 600	74 - 79
<b>LT</b>	Naudojimo instrukcijos – SARSTEDT „Multivette® 600“	80 - 85
<b>LV</b>	Lietošanas norādes – SARSTEDT Multivette® 600	86 - 91
<b>NL</b>	Gebruiksaanwijzing – SARSTEDT Multivette® 600	92 - 97
<b>NO</b>	Bruksanvisning – SARSTEDT Multivette® 600	98 - 103
<b>PL</b>	Instrukcja obsługi – Multivette® 600 SARSTEDT	104 - 109
<b>PT</b>	Instruções de utilização – Multivette® 600 SARSTEDT	110 - 115
<b>RO</b>	Instrucțiuni de utilizare – Multivette® 600 SARSTEDT	116 - 121
<b>RU</b>	Инструкция по применению – SARSTEDT Multivette® 600	122 - 127
<b>SK</b>	Návod na Použitie – SARSTEDT Multivette® 600	128 - 133
<b>SL</b>	Navodila za uporabo – Multivette® 600 SARSTEDT	134 - 139
<b>SV</b>	Bruksanvisning – SARSTEDT Multivette® 600	140 - 145
<b>TH</b>	คำแนะนำในการใช้งาน – ระเบียบการจะเล็ด Multivette® 600 ของ SARSTEDT	146 - 151
<b>TR</b>	Kullanım için talimatlar – SARSTEDT Multivette® 600	152 - 157
<b>ZH</b>	使用说明 – SARSTEDT Multivette® 600	158 - 163

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## Intended use

The Multivette® 600 is a blood collection system and is used for venous blood collection with Luer needles or for capillary blood collection using the end-to-end principle. The Multivette® 600 is used for collecting, transporting and processing blood samples in the clinical laboratory. The product is intended for use in a professional environment by qualified medical and laboratory personnel.

## Product description

The Multivette® 600 consists of a plastic tube, a colour-coded plastic screw cap and various additives (preparations). The selection of the additive depends on the analytical test method. It is specified by the manufacturer of the test reagents and/or the analysis equipment on which the test is conducted.

Multivette® 600 screw cap colour codes:

Additive	Letter code	Cap colour based on ISO 6710*	Cap colour based on BS 4851**
<b>Serum tube</b>			
Multivette® 600 with serum	CAT	red	white
Multivette® 600 with serum gel	CAT	brown	brown
<b>Lithium heparin tube</b>			
Multivette® 600 with lithium heparin	LH	green	orange
Multivette® 600 with lithium heparin gel	LH	green	orange
<b>EDTA tube</b>			
Multivette® 600 with K3 EDTA	K3E	violet	red
<b>Tube with glycolysis inhibitors</b>			
Multivette® 600 with fluoride/heparin	FH	-	yellow

\*DIN EN ISO 6710: Single-use containers for human venous blood specimen collection.

\*\*British Standard BS4851:1982 Specification for single use labelled medical specimen containers for haematology and biochemistry, revoked.

## SARSTEDT Multivette® 600 Serum / Serum Gel

The Multivette® 600 with serum contains a coagulation activator and is used for collecting serum. The Multivette® 600 with serum gel contains a coagulation activator and is used for collecting serum. After centrifugation, the gel forms a layer separating the coagulum and the serum. Multivettes 600 with serum / serum gel are intended for clinical chemistry tests on serum.

## SARSTEDT Multivette® 600 Lithium Heparin

The Multivette® 600 with lithium heparin is coated with lithium heparin as an anticoagulant and is used for collecting plasma. The Multivette® 600 with lithium heparin gel contains lithium heparin as an anticoagulant and is used for collecting plasma. After centrifugation, the gel forms a layer separating the plasma and the cellular blood components. The Multivettes 600 with lithium heparin / lithium heparin gel are intended for clinical chemistry tests on plasma.

## SARSTEDT Multivette® 600 K3 EDTA

The Multivette® 600 with K3 EDTA is coated with K3 EDTA as an anticoagulant and is intended for haematological tests.

## SARSTEDT Multivette® 600 Fluoride/Heparin

The Multivette® 600 with fluoride/heparin is coated with the anticoagulant heparin and contains fluoride for glycolysis inhibition. It is used to test glucose levels.

**NOTE:** *Fluoride (used for measuring glucose) can cause an increase in haemolysis. Further information on substances that can cause interference can be found in the respective instructions for use of the assay manufacturer.*

## Safety information and warnings

1. General precautions: Use gloves and general personal protective equipment to protect yourself from possible exposure to potentially infectious specimen material and transmitted pathogens.
2. Handle all biological specimens and collection utensils according to the guidelines and procedures in your facility. Consult a doctor in the event of direct contact with biological samples, as they can transmit infectious diseases.
3. Dispose of all sharp objects (needles) for blood collection in appropriate disposal containers.
4. Under- or overfilling the Multivette® 600 leads to an incorrect ratio of blood to preparation/additive and can result in false analysis results.
5. The blood collected and processed with the Multivette® 600 is not intended to be re-injected into the human body.
6. Do not use the products after the expiry date. The expiry date of the product is the last day of the month and year as indicated.

## Storage

Store the product at room temperature.

## Limitations

1. If blood samples are stored in the Multivette® 600, the stability of the analytes should be evaluated by the respective laboratory or taken from the instructions for use of the analysis device manufacturer.
2. Even if plasma or serum is separated by centrifuging the Multivette® 600 and/or a barrier is present, all cells are not necessarily completely separated. Residual metabolism or natural degradation can affect the concentrations of analytes. The stability of the analyte should be evaluated with respect to the storage containers and the conditions of the respective laboratory.
3. For therapeutic medication, the suitability of the sample material must be checked in the manufacturer's instructions for use for the assay / analysis device.

## Specimen collection and handling

**READ ALL OF THIS DOCUMENT BEFORE STARTING THE BLOOD COLLECTION TRANSFER.**

## Materials needed for collecting a specimen/blood

### Material required but not supplied for the blood collection:

1. All necessary Multivettes 600, labelled by size and additive.
2. Gloves, gown, goggles or other suitable protective wear to protect from bloodborne pathogens or potentially infectious materials.
3. Labels for specimen identification.
4. (Safety) lancet for capillary blood collection and (safety) needles or Safety-Multifly® Needles for venous blood collection.
5. Disinfectant to clean the collection site (follow guidelines of the facility for sterile specimen collection to prepare the collection site). Do not use any alcohol-based cleaning materials if the specimens are to be used for a blood alcohol test.
6. Dry, clean swabs.
7. Plasters.
8. Sharps disposal container for the safe disposal of used material.

## Recommended order of collection for venous blood collection (CLSI\* GP39):

1. Tubes for coagulation tests of serum with/without gel
2. Tubes with heparin with/without gel
3. Tubes with EDTA with/without gel
4. Tubes with glycolysis inhibitors
5. Tubes with other additives

**NOTE:** Follow the regulations of your institution for the order of collection.

## Recommended order of collection for capillary blood collection (CLSI\* GP42):

1. Tubes with EDTA
2. Tubes with lithium heparin / lithium heparin gel
3. Tubes with fluoride
4. Tubes for coagulation tests serum / serum gel


**NOTE:** The first drop of blood is discarded.

**NOTE:** Follow the regulations of your institution for the order of collection.

## Venous blood collection with a Luer needle


### General instructions:

WEAR GLOVES TO MINIMISE THE RISK OF EXPOSURE DURING BLOOD COLLECTION.


- 1 For simple venous collections, attach a Luer needle to the capillary tube. The special design enables automatic filling of the inner tube using only the venous pressure.


Stop collecting blood once the blood flow has reached the end of the capillary tube (600 µl).


- 2 After the venous blood collection is complete, first remove the needle


- 3 ... and dispose of it.

- 4 Immediately after that, drain the capillary tube by pulling it out and dispose of it as a single unit.

- 5 Seal the container using the screw cap provided...

- 6 ... and mix the specimen carefully by inverting it gently.

### Information on coagulation / on the Multivette® 600 Serum:


After collecting blood, let the blood coagulate in the Multivette® 600 Serum / Serum Gel for 30 minutes in an upright position to ensure a clear layer of separation after centrifugation.

The recommended times are based on an intact coagulation process. Blood from patients with coagulation disorders due to their illness or under anticoagulant therapy requires more time for complete coagulation.



## Capillary technique using the end-to-end principle

### General instructions:

WEAR GLOVES TO MINIMISE THE RISK OF EXPOSURE DURING CAPILLARY BLOOD COLLECTION.


- 1 Hold the Microvette® 600 horizontally or tilted slightly and collect the drops of blood with the capillary tube. The blood collection is complete once the blood flow has reached the end of the capillary tube (600 µl).

- 2 Immediately after that, drain the capillary tube by pulling it out and dispose of it as a single unit.

- 3 Seal the container using the screw cap provided...

- 4 ... and mix the specimen carefully by inverting it gently.

#### Information on coagulation / on the Multivette® 600 Serum:

After collecting blood, let the blood coagulate in the Multivette® 600 Serum / Serum Gel for 30 minutes in an upright position to ensure a clear layer of separation after centrifugation.

The recommended times are based on an intact coagulation process. Blood from patients with coagulation disorders due to their illness or under anticoagulant therapy requires more time for complete coagulation.

## Centrifugation

Refer to the centrifuge instructions for use for information on disinfecting the centrifuge.

### IMPORTANT!

Centrifuging Multivettes 600 with cracks or centrifuging at an excessive centrifugal force can lead to breakage of the Multivettes 600, which can release potentially hazardous substances. Please see our centrifugation recommendations for the various additives in the table below.

Centrifuge inserts are to be selected according to the size of the Multivettes 600 used. The relative centrifugal force has the following ratio to the selected rpm:

$$\text{RCF (g-force)} = 11.2 \times r \times (\text{rpm}/1,000)^2,$$

"RCF": "Relative centrifugal force",





"rpm": "Revolutions per minute" (min<sup>-1</sup>),

"r" [in cm]: "Radius of rotation from the middle of the centrifuge to the bottom" of the Multivette® 600.

Multivettes 600 without gel can be centrifuged in centrifuges with a fixed-angle or swing-out rotor. Multivettes 600 with gel are designed solely for use in centrifuges with swing-out rotor. Centrifugation in centrifuges with a fixed-angle rotor is not validated by SARSTEDT and is not recommended. In addition, centrifugation should be conducted no later than 2 hours after collection, as this could otherwise yield false results.

Multivettes 600 should be centrifuged in accordance with the centrifugation conditions listed below. If other conditions are used, they must be validated by the user. It must be ensured that the Multivettes® 600 fit well in the centrifuge inserts. Multivettes 600 that protrude above the insert can get caught in the centrifuge head and be broken. Ensure that the centrifuge is filled evenly. Please refer to the centrifuge instructions for use for further information.

**CAUTION! Do not remove broken Multivettes 600 by hand.**

Image	Additive	Times and g values (=RCF*)			
		Standard recommendation		Alternative recommendation	
	Serum	5 min	10,000 x g	or	
	Lithium heparin	5 min	2,000 x g	or	10 min 2,000 - 10,000 x g
	Glucose fluoride	5 min		or	
	EDTA <sup>2</sup>	5 min	10,000 x g	or	
	Serum gel <sup>1</sup>	5 min	10,000 x g	or	10 min 4,000 - 10,000 x g
	Lithium heparin gel <sup>1</sup>				

<sup>1</sup> For tubes prepared with gel, we recommend using swing-out rotors.

<sup>2</sup> Only for plasma separation, not for haematological tests.

\*RCF: Relative centrifugal force

Centrifugation at 20°C

## Disposal

1. The general hygiene guidelines and regulations for the proper disposal of infectious material must be observed and complied with.
2. Disposable gloves prevent the risk of infection.
3. Contaminated or filled blood collection systems must be disposed of in suitable containers for hazardous biological waste that can be subsequently autoclaved and incinerated.
4. They must be disposed of in a suitable incinerator or by autoclaving (steam sterilisation).

## Product-specific standards and guidelines as amended/updated

### Product-specific standards:

CLSI\* GP39 "Tubes and Additives for Venous and Capillary Blood Specimen Collection", Approved Standard.

CLSI\* GP41 "Collection of Diagnostic Venous Blood Specimens", Approved Standard.

CLSI\* GP42 "Collection of Capillary Blood Specimens", Approved Standard.

CLSI\* GP44 "Procedures for the Handling and Processing of Blood Specimens for Common Laboratory Tests", Approved Guideline.

WHO/DIL/LAB/99.1 Rev02 "WORLD HEALTH ORGANIZATION, et al. Use of anticoagulants in diagnostic laboratory investigations. Geneva: World Health Organization, 2002".

\*CLSI (Clinical and Laboratory Standards Institute)

### Key for symbols and labels:

	Article number
	Batch number
	Use by
	CE marking
	<i>In vitro</i> diagnostic device
	Follow the instructions for use
	If reused: Risk of contamination
	Store away from sunlight
	Store in a dry place
	Manufacturer
	Country of manufacture

Technical modifications reserved.

All serious incidents relating to the product must be reported to the manufacturer and the competent national authority.

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\*\*British Standard BS4851:1982 Specification for single use labelled medical specimen containers for haematology and biochemistry, отменен.

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- CLSI\* GP39 "Tubes and Additives for Venous and Capillary Blood Specimen Collection", Approved Standard.
  - CLSI\* GP41 "Collection of Diagnostic Venous Blood Specimens", Approved Standard.
  - CLSI\* GP42 "Collection of Capillary Blood Specimens", Approved Standard.
  - CLSI\* GP44 "Procedures for the Handling and Processing of Blood Specimens for Common Laboratory Tests", Approved Guideline.
  - WHO/DIL/LAB/99.1 Rev02 "WORLD HEALTH ORGANIZATION, et al. Use of anticoagulants in diagnostic laboratory investigations. Geneva: World Health Organization, 2002".
- \*CLSI (Clinical and Laboratory Standards Institute)

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\*\*British Standard BS4851:1982 Specification for single use labelled medical specimen containers for haematology and biochemistry, atšauktas.

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CLSI\* GP39 " Tubes and Additives for Venous and Capillary Blood Specimen Collection" , Approved Standard.

CLSI\* GP41 " Collection of Diagnostic Venous Blood Specimens" , Approved Standard.

CLSI\* GP42 " Collection of Capillary Blood Specimens" , Approved Standard.

CLSI\* GP44 " Procedures for the Handling and Processing of Blood Specimens for Common Laboratory Tests" , Approved Guideline.

WHO/DIL/LAB/99.1 Rev02 " WORLD HEALTH ORGANIZATION, et al. Use of anticoagulants in diagnostic laboratory investigations. Geneva: World Health Organization, 2002" .

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