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SARSTEDT Stuhlprobenröhre

DE	Gebrauchshinweis – SARSTEDT Stuhlprobenröhre	2 - 6
EN	Instructions for Use – SARSTEDT faeces tube	7 - 11
BG	Инструкции за употреба – SARSTEDT епруветка за фекални проби	12 - 16
CS	Návod k obsluze – Zkumavka na vzorek stolice SARSTEDT	17 - 21
DA	Brugsanvisning – SARSTEDT fæcesprøverør	22 - 26
EL	Οδηγίες χρήσεως – Σωληνάριο δείγματος κοπράνων SARSTEDT	27 - 31
ES	Instrucciones de uso – Tubo para heces SARSTEDT	32 - 36
ET	Kasutusjuhend – Roojaproovikatsutid SARSTEDT	37 - 41
FR	Mode d'emploi – Tube de recueil de selles SARSTEDT	42 - 46
HR	Uputa za upotrebu – SARSTEDT epruveta za uzimanje uzoraka stolice	47 - 51
HU	Használati utasítás – SARSTEDT székletminta-cső	52 - 56
IT	Istruzioni d'uso – Provetta per campioni di feci SARSTEDT	57 - 61
KO	사용 설명서 – SARSTEDT 대변 검체 튜브	62 - 66
LT	Naudojimo instrukcijos – SARSTEDT išmatų mėginių mėgintuvėliai	67 - 71
LV	Lietošanas norādes – SARSTEDT fēču paraugu caurulīte	72 - 76
NL	Gebruiksaanwijzing – SARSTEDT-buisjes voor faecesmonsters	77 - 81
NO	Bruksanvisning – SARSTEDT avføringsprøverør	82 - 86
PL	Instrukcja obsługi – próbówka na kał SARSTEDT	87 - 91
PT	Instruções de utilização – Tubo de amostra para colheita de fezes SARSTEDT	92 - 96
RO	Instrucțiuni de utilizare – Tuburi de coprocultură SARSTEDT	97 - 101
RU	Инструкция по применению – Контейнеры для сбора кала SARSTEDT	102 - 106
SK	Návod na Použitie – SARSTEDT skúmavka na vzorky stolice	107 - 111
SL	Navodila za uporabo – posodica za vzorec blata SARSTEDT	112 - 116
SV	Bruksanvisning – SARSTEDT Faecesrör	117 - 121
TH	คำแนะนำในการใช้งาน – หลอดเก็บตัวอย่างอุจจาระของ SARSTEDT	122 - 126
TR	Kullanım için talimatlar – SARSTEDT gaita kabı	127 - 131
ZH	使用说明 – SARSTEDT 粪便样本管	132 - 136

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Intended use


Faeces tubes are sample containers used for the simple and safe collection, storage and processing of a stool sample for *in vitro* diagnostic determinations. They are available in different sizes and are closed with either a screw cap or a push cap. A spoon is mounted in the lids, with which the stool sample can be picked up simply and cleanly and transferred into the tube.

The products are intended for sample collection under the instruction of medical professionals and for use and application in a professional environment by medical professionals and laboratory personnel.

Product description

The faeces tubes consist of a clear plastic container with a matching screw cap or push cap and are available in different filling volumes with different opening diameters. Depending on the design, the faeces tubes are made of different materials and differ in dimensions, stool spoon type and bottom shape. Depending on the version, they are available with a label, imprint or suitable shipping container. Sterile and non-sterile versions are available.

Product overview

Faeces tubes:

Type	Description
A	Faeces tube with screw cap, length/ø 107/25 mm
B	Faeces tube with screw cap, length/ø 101/16.5 mm
C	Faeces tube with screw cap, length/ø 76/20 mm
D	Faeces tube with push cap, length/ø 75/23.5 mm
E	Faeces tube with screw cap for obtaining a defined stool sample (approx. 1 ml), with round bottom with skirted base, length/ø 101/16.5 mm, with enclosed spatula.

Safety information and warnings

1. Check each individual package for damage before use and do not use if packaging is damaged.
2. Follow the directions for handling.
3. General precautions: Use gloves and general personal protective equipment to protect yourself from possible exposure to potentially infectious specimen material and transmitted pathogens.
4. Handle all biological specimens and collection utensils according to the guidelines and procedures in your facility. In the event of contact with biological samples, consult a doctor, as this can transmit infectious diseases. Always follow the safety guidelines and procedures of your facility.
5. Always store the closed, filled faeces tubes in an upright position, not lying down.
6. Faeces tubes are intended for single use. Dispose of all products in the hazardous biological waste containers.
7. Do not use the faeces tubes after the expiry date. The expiry date of the faeces tubes is the last day of the month and year indicated.

Storage

Store the products at room temperature.

Transport

Version B faeces tubes with screw cap, together with an ADR/IATA (packaging instruction P650) compliant secondary container, incl. suction insert and type-tested outer packaging, comply with the ADR P650/IATA regulations.

Always transport faecal samples to the laboratory as quickly as possible.

Limitations

1. The storage period and temperature are dependent on the expiry date of the parameters for investigation. Clarify the time within which the stool analysis should be carried out with the responsible laboratory or refer to the analysis manufacturer's instructions for use.
To potentially prolong sample stability for some analytes, you can cool the sample material at 2–8 °C.
2. Whenever possible, take samples before initiating pathogen-damaging measures (antibiotics, antiseptics, antiviral therapeutics).
In the case of pathogen-damaging measures, clarify the suitability of the sample material with the competent laboratory.
3. Homogenise the stool sample by gentle stirring prior to examination.

Collecting and handling specimens

READ THIS DOCUMENT COMPLETELY BEFORE STARTING WITH THE STOOL SAMPLE COLLECTION.

Equipment required for specimen collection

1. Faeces tube
2. If necessary, Stool catcher
3. Gloves, gowns, eye protection or other suitable protective clothing to protect against pathogens or potentially infectious materials transmitted by sample material.
4. If necessary, Label
5. If necessary, suitable ADR/IATA-compliant shipping containers and type-tested outer packaging

Performing the stool sample collection

General information

In some cases, 2 to 3 separate stool samples will need to be sent for examination.

This increases the detection rate in the case of discontinuous pathogen excretion - a single negative microbiological test result does not exclude an infection!

If several samples are collected at specific intervals, the samples must be individually forwarded to the laboratory immediately, not kept for sending later!

In a person has several bowel movements per day, the sample should, if possible, be collected during the first bowel movement.



When obtaining a stool sample, it is important that there be no contact with toilet flushing water and other contaminants. Only in this way can a meaningful test result be ensured.

If you have a washout toilet, remove the water before depositing the stool.

Alternatively, you can – especially in the case of a washdown toilet – stretch a so-called "stool catcher" over the toilet seat.

A stool catcher avoids:

- Dilution of the stool sample by flushing water or urine
- Contact of the stool sample with toilet cleaner or fragrances

The different versions of the faeces tubes are described below:

A: When using a faeces tube with a screw cap


B: When using a faeces tube with push cap

C: When using a faeces tube with a screw cap to obtain a defined stool sample (approx. 1 ml)

A: When using a faeces tube with a screw cap

After the patient has deposited the stool sample:


1. Open the lid of the faeces tube by turning it anticlockwise (fig. 1).
2. Remove material from 3 to 5 different parts of a stool sample using the stool spoon integrated into the lid. Please collect in particular from locations containing mucus and blood (fig. 2a & 2b). The amount of stool required may vary depending on the examination. Follow the instructions of the responsible analysis laboratory; if applicable, the label includes a note.
3. Close the faeces tube tightly with the screw cap by turning it clockwise. The spoon remains in the container, please DO NOT remove (fig. 3).
4. Mark the faeces tube with the required patient data (e.g. name, date of birth and date of collection). Use the already applied label for this purpose, or, if using a faeces tube without a label, a separate label must be applied to the tube.
5. Send each stool sample to the laboratory as soon as possible and store it in a cool place until delivered, but do not freeze. (fig. 4)
6. Store in a cool place in the laboratory until processing in order to avoid overgrowth of the enteropathogenic pathogens by the physiological faecal flora.

B: When using a faeces tube with push cap

After the patient has deposited the stool sample:


1. Open the cap of the faeces tube by pulling it out upwards (fig. 1).
2. Remove material from 3 to 5 different locations of a stool sample using the stool spoon integrated into the cap. Please collect in particular from locations containing mucus and blood (fig. 2a & 2b). The amount of stool required may vary depending on the examination. Follow the instructions of the responsible analysis laboratory; if applicable, the label includes a note.
3. Close the faeces tube with the cap by pressing it tightly into the tube. The spoon remains in the container, please DO NOT remove (fig. 3).
4. Mark the faeces tube with the required patient data (e.g. name, date of birth and date of collection). Use the already applied label for this purpose, or, if using a faeces tube without a label, a separate label must be applied to the tube.
5. Send each stool sample to the laboratory as soon as possible and store it in a cool place until delivered, but do not freeze (fig. 4).
6. Store in a cool place in the laboratory until processing in order to avoid overgrowth of the enteropathogenic pathogens by the physiological faecal flora.

C: When using a faeces tube with a screw cap to obtain a defined stool sample (approx. 1 ml)

After the patient has deposited the stool sample:

1. Open the lid of the faeces tube by turning it anticlockwise (fig. 1).
2. Remove material from 3 to 5 different parts of a stool sample using the stool spoon integrated into the lid. Please collect in particular from locations containing mucus and blood. Make sure to fill the entire spoon (fig. 2a & 2b).
3. Scrape off an excess with the enclosed spatula (fig. 3).
4. Discard the excess stool and spatula. Please observe your in-house instructions and the disposal instructions below.
5. Close the faeces tube tightly with the screw cap by turning it clockwise. The spoon remains in the container, please DO NOT remove (fig. 4).
6. Mark the faeces tube with the required patient data (e.g. name, date of birth and date of collection). Use the already applied label for this purpose, or, if using a faeces tube without a label, a separate label must be applied to the tube.
7. Send each stool sample to the laboratory as soon as possible and store it in a cool place until delivered, but do not freeze (fig. 5).
8. Store in a cool place in the laboratory until processing in order to avoid overgrowth of the enteropathogenic pathogens by the physiological faecal flora.
9. To prepare a suspension, open the lid of the faeces tube by turning it anticlockwise (fig. 6).
10. Add 2 ml of distilled water (fig. 7).
11. Close the faeces tube tightly with the screw cap by turning it clockwise. The spoon remains in the container, please DO NOT remove (fig. 8).
12. Mix the faeces tube simply by shaking it (fig. 9).
13. Centrifuge the faeces tube to produce a supernatant.
14. The supernatant can now be sent for analysis.

Centrifugation

SARSTEDT faeces tube version E is designed for a maximum of 3,500 × g for 15 minutes (RCF: relative centrifugal force – g-force). The relative centrifugal force has the following ratio to the selected rpm:

$$\text{RCF (g-force)} = 11.2 \times r \times (\text{rpm}/1000)^2,$$

"RCF": "Relative centrifugal force", RCF

"rpm": "Revolutions per minute",

"r" [in cm]: "Centrifuge radius from the centre of the centrifuge to the bottom of the faeces tube.

Only suitable centrifuge inserts should be used.

IMPORTANT!

Centrifugation of faeces tubes with cracks or centrifugation if the centrifugal acceleration is too high can lead to breakage of the faeces tubes, potentially releasing hazardous substances. Faeces tubes should be centrifuged according to the centrifugation conditions listed above. If other conditions are used, they must be validated by the respective user.

Ensure that the faeces tubes fit into the centrifuge inserts. Faeces tubes that protrude beyond the insert can catch on the centrifuge head and be broken. Ensure that the centrifuge is filled evenly. Refer to the instructions for use of the centrifuge for this purpose.

CAUTION! Do not remove broken faeces tubes by hand.

Refer to the centrifuge operating manual for instructions on disinfecting the centrifuge.

Disposal

1. The general hygiene guidelines and regulations for the proper disposal of infectious material must be observed and complied with.
2. Disposable gloves prevent the risk of infection.
3. Contaminated or filled faeces tubes and their accessories must be disposed of in suitable containers for hazardous biological waste that can be subsequently autoclaved and incinerated.
4. They must be disposed of in a suitable incinerator or by autoclaving (steam sterilisation).

Key for symbols and labels:

Article number

Batch number

Use by

CE marking

In vitro diagnostic device

Follow the instructions for use

If reused: risk of contamination

Keep away from sunlight

Store in a dry place

Manufacturer

Country of manufacture

The following also applies to sterile products:

Product sterile inside

Do not resterilise

Single sterile barrier system with external protective packaging

Technical modifications reserved.

All serious incidents relating to the product must be reported to the manufacturer and the competent national authority.

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