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Intended use

Prepared sample tubes for blood glucose testing from haemolysate are used with End-to-End capillaries and a haemolysis solution. The End-to-End capillary is used for capillary blood collection. Prepared sample tubes are used for holding an End-to-End capillary filled with capillary blood and for the haemolysis and processing of blood samples in clinical laboratories. The products are intended for use in a professional environment by qualified medical and laboratory personnel.

Product description

Prepared sample tubes for blood glucose testing from haemolysate consist of a plastic tube with a conical false bottom and a flat bottom, a white plastic screw cap or a grey plastic membrane screw cap, various additives (preparations) and a separate plastic capillary with various additives (preparations).

Colour code	Pre dosing	Dosing quantity/ capillary	Capillary prep.	Length / Ø in mm
White	Haemolysis solution for blood glucose testing (SARSTEDT)	1,000 µl/20 µl	K2-EDTA	92 / 15.3 round bottom
Grey	Haemolysis solution for blood glucose testing (SARSTEDT)	2,000 µl/10 µl	Ammonium heparin	66 / 11.5 round bottom

Safety information and warnings

1. General precautions: Use gloves and other general personal protective equipment to protect yourself against blood and possible exposure to pathogens transmitted by blood.
2. Handle all biological specimens and sharp blood collection utensils (needles) according to the guidelines and procedures of your facility. In case of direct contact with biological specimens or of a needlestick injury, consult a doctor due to the risk of transmission of HIV, HCV, or HBV or other infectious diseases. Always follow the safety guidelines and procedures of your facility.
3. Dispose of all sharp objects (needles) for blood collection in suitable sharps containers.
4. under- or overfilling the prepared sample tubes leads to an incorrect ratio of blood to preparation/additive and can give rise to false analysis results.
5. Do not use the products after the expiry date. The expiry date of the product is the last day of the month and year indicated.

Storage

Store the products at room temperature.

Transport

The products are primary receptacles according to the ADR (Packing instruction P650) and IATA regulations.

Limitations

1. The storage period and temperature of a filled sample tube depend on the stability of the analytes for investigation. The assessment is made by the laboratory, or information is obtained from the manufacturer's instructions for use for the test reagents / analysis device.
2. For therapeutic medication, the suitability of the sample material must be checked in the manufacturer's instructions for use for the test reagents / analysis device.

Collecting and handling specimens

READ THIS DOCUMENT COMPLETELY BEFORE STARTING TO DRAW BLOOD.


Materials needed for collecting a specimen/blood:

1. The required prepared tubes and the corresponding End-to-End capillaries.
2. Gloves, gown, goggles or other suitable protective gear to protect from bloodborne pathogens or potentially infectious materials.
3. Labels for specimen identification.
4. (Safety) lancet or incision lancet.
5. Disinfectant to clean the collection site (follow guidelines of the facility for sterile specimen collection to prepare the collection site). Do not use any alcohol-based cleaning materials if the specimens are to be used for a blood alcohol test.
6. Dry, clean swabs.
7. Plasters.
8. Sharps disposal container for the safe disposal of used material.

Capillary technique with the End-to-End capillary

General instructions:

WEAR GLOVES TO MINIMISE THE RISK OF EXPOSURE DURING CAPILLARY BLOOD COLLECTION.

- 1 Prepared sample tubes for blood glucose testing from haemolysate consisting of a prepared sample tube and End-to-End capillary.
- 2 Clean the puncture site with disinfectant. Puncture the skin in a way that ensures good blood flow. Discard the first drop. Then hold the End-to-End capillary horizontally with one end in the middle of the drop of blood and completely fill the End-to-End capillary without air bubbles.
- 3 Open the prepared sample tube and place the End-to-End capillary filled with blood in it.
- 4 Close the prepared sample tube and swirl it gently.

Note:

For the stability of the analytes, follow the instructions for use of the test reagents / analysis device manufacturer.

Disposal

1. The general hygiene guidelines and regulations for the proper disposal of infectious material must be observed and complied with.
2. Disposable gloves prevent the risk of infection.
3. Contaminated or filled prepared sample tubes must be disposed of in suitable containers for hazardous biological waste that can be subsequently autoclaved and incinerated.
4. They must be disposed of in a suitable incinerator or by autoclaving (steam sterilisation).

Product-specific standards and guidelines as amended/updated

CLSI* GP44 "Procedures for the Handling and Processing of Blood Specimens for Common Laboratory Tests", Approved Guideline .
 WHO/DIL/LAB/99.1 Rev02 "WORLD HEALTH ORGANIZATION, et al. Use of anticoagulants in diagnostic laboratory investigations. Geneva: World Health Organization, 2002".
 *CLSI (Clinical and Laboratory Standards Institute)

Key for symbols and labels:

	Article number
	Batch number
	Use by
	CE marking
	<i>In vitro</i> diagnostic device
	Follow the instructions for use
	If reused: Risk of contamination
	Store away from sunlight
	Store in a dry place
	Manufacturer
	Country of manufacture

Technical modifications reserved.

All serious incidents relating to the product shall be notified to the manufacturer and the competent national authority.

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