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SARSTEDT Sediplus® S 2000 NX

DE	Gebrauchshinweis – SARSTEDT Sediplus® S 2000 NX	2 - 55
EN	Instructions for Use – SARSTEDT Sediplus® S 2000 NX	56 - 109
CS	Návod k obsluze – SARSTEDT Sediplus® S 2000 NX	110 - 163
ES	Instrucciones de uso – SARSTEDT Sediplus® S 2000 NX	164 - 217
FR	Mode d'emploi – SARSTEDT Sediplus® S 2000 NX	218 - 271
IT	Istruzioni d'uso - SARSTEDT Sediplus® S 2000 NX	272 - 325
NL	Gebruiksaanwijzing – SARSTEDT Sediplus® S 2000 NX	326 - 379
PL	Instrukcja obsługi – SARSTEDT Sediplus® S 2000 NX	380 - 433
PT	Instruções de utilização – Sediplus® S 2000 NX SARSTEDT	434 - 487
SK	Návod na použitie – SARSTEDT Sediplus® S 2000 NX	488 - 541

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Basic notes!

Before initial operation of the Sediplus® S 2000 NX, observe these Instructions for Use!
 The basic prerequisite for the correct handling and trouble-free operation of the device is knowledge of the contents of these Instructions for Use.

Online instructions for use: Sediplus® S 2000 NX; GB 557 at:
www.sarstedt.com/en/download/instructions-for-use

Copyright:

SARSTEDT AG & Co. KG is the copyright holder of these Instructions for Use.
 The Instructions for Use are intended only for the operating personnel and for the purchaser of the device. These Instructions for Use may not be reproduced or distributed in whole or in part without the written consent of SARSTEDT AG & Co. KG. Violations can have criminal consequences.

Please keep the Instructions for Use as a reference for information on your device.

Technical modifications reserved!

Sediplus®, Sedivette®, Sarmix® and Monovette® are registered trademarks of SARSTEDT AG & Co. KG

Nümbrecht, May 2022
 SARSTEDT AG & Co. KG

<p>Manufacturer and customer service address:</p>	<p>Device data: (to be completed by the customer)</p>
 <p>SARSTEDT AG & Co. KG PO Box 1220 D-51502 Nümbrecht</p> <p>Phone: +49 (0) 22 93-30 50 Fax: +49 (0) 22 93-305 282 E-Mail: info@sarstedt.com www.sarstedt.com</p>	<p>Type: Sediplus® S 2000 NX</p> <p>Serial No.: Place of installation: Issue date: Inventory No.: Software version (control panel): Software version (measuring table electronics):</p>

Last modified: 08/08/2022
 Valid from serial number: BSG0001

Table of Contents

Basic notes!	56
1. General information	59
1.1 Safety information	59
1.2 Intended use	63
1.3 EC Declaration of Conformity	63
2. Introduction	64
2.1 Description	64
2.2 The "blood sedimentation" measurement	64
2.3 Setting up the device / Installation	66
2.3.1 Technical data:	67
2.3.2 Unpacking, transport and storage	69
2.3.3 Delivery contents.....	69
2.3.4 Connecting a PC-AT keyboard and / or a barcode scanner.....	70
2.3.5 Connecting an optional printer.....	71
2.3.6 Description of the HL7 interface	72
3. The functional elements of the Sediplus® S 2000 NX	74
3.1 The measuring plate.....	74
3.2 The light emitting diodes of the measuring plate	75
3.3 The Display	76
4. Operation of the Sediplus® S 2000 NX	77
4.1 Switching the device on	77
4.2 Preparation of samples.....	78
4.3 Inserting the S-Sedivette® and starting a measurement	79
4.4 Selecting the measurement times.....	82
4.5 Entering the ID number.....	82
4.5.1 Entering the ID number using a barcode scanner (optional)	83
4.5.2 Entering the ID number via the touch display	83
4.6 Viewing already determined values during a measurement	84
4.7 Acoustic support.....	85
4.8 Reprinting / repeating the printing process	86
4.9 Deleting measurement results	86
5. Menu	87
5.1 Settings.....	87
5.1.1 Language.....	87
5.1.2 Start ID.....	88
5.1.3 Measurement settings	88
5.1.3.1 Measurement type	88
5.1.3.2 "Underfill limit 5 ... 10 mm"	89
5.1.4 Deleting measurement results	89
5.1.5 Date	90
5.1.6 Time.....	90
5.1.7 Printer settings	91
5.1.8 Network settings / HL7.....	92
5.2 Service.....	93
5.3 Errors	93

6.	Faults and operating errors.....	93
6.1	Fault in the opening and closing of the measuring plate.....	93
6.2	Underfilled S-Sedivette®	94
6.3	Incorrect measurement	95
6.4	Premature removal of an S-Sedivette®	96
6.5	Mains failure	98
7.	Service and maintenance	98
7.1	Cleaning the device	98
7.2	Quality controls with test Sedivettes	99
7.3	Measurement optics.....	101
7.4	Maintenance at the manufacturer's premises	101
8.	Decommissioning / Disposal	102
9.	General precautions:.....	102
10.	Limitations	102
11.	Order information	103
12.	Warranty and liability.....	104
13.	Glossary.....	104
14.	List of illustrations	105
15.	Literature.....	106
16.	Key for symbols and labels:.....	107
17.	Contamination questionnaire	109

1. General information

1.1 Safety information

The Sediplus® S 2000 NX blood sedimentation measuring device has been built according to the latest engineering standards and the recognised safety regulations. Take note of the relevant safety information and guidelines as well as occupational safety and accident prevention regulations for use in the laboratory.

Nevertheless, risks for users and third parties as well as damage to the device or other material assets can still arise during operation.

The device should be used only:

- For the intended purpose.
- In faultless condition in terms of safety.

During the device's warranty period, only the company SARSTEDT AG & Co. KG or persons authorised by the company SARSTEDT AG & Co. KG may repair the device.

Improper handling or repairs will render any warranty claims void.

- Measurement results must always be checked for plausibility by the operating personnel, the laboratory management or the attending doctor.
- The general regulatory guidelines and statutory regulations for properly determining blood sedimentation must be noted and complied with.

The device operates at a voltage of 110-230 V, 50-60 Hz. No water or other liquids may enter the device.

- When connecting the device to the mains power, you may use only the supplied power cords.
- Check to see whether your mains voltage matches the information on the type label. The device may only be connected to an earthed shockproof socket.
- Check for damage to the connection cable before start-up. Damaged cables may not be used under any circumstances.
- Check the device for external damage before start-up. Never work with a damaged device.
- Remove the plug from the socket if the device will not be used for a long period of time or if it needs to be cleaned.

The S-Sedivettes and the blood samples in the S-Sedivettes must be handled and disposed of in the proper manner. Note the applicable regulations in your country and follow the safety provisions in your laboratory.

The device may be operated or complemented only using the accessories described in this manual.

OBLIGATION OF THE OPERATOR

The device's operator undertakes to commission only persons who have read and understood this manual with working on the device. This should be confirmed with their signature.

Checks to ensure personnel are aware of workplace safety should be carried out at regular intervals.

OBLIGATION OF THE PERSONNEL

Persons working with the device undertake to:

- Read through these Instructions for Use before commencing work.
- Comply with the applicable accident prevention regulations and requirements.

Both should be confirmed by the recipient's signature.

As there are many error sources that can affect the Erythrocyte Sedimentation Rate (ESR), it is imperative that the information below be noted when setting up the device:

- Vibrations of the filled sedimentation system lead to increased sedimentation values. The device must therefore be placed on a fixed table or on a vibration-free surface. For example, the device must not stand on a laboratory table together with a centrifuge.
- The S-Sedivettes may not be scratched, dusty or contaminated under any circumstances.
- Setting up the filled lowering system on an incline leads to increased sedimentation values. The device must therefore be set up in such a way that the S-Sedivettes can be loaded from above and the S-Sedivettes stand upright.
- The measuring plate will move upwards during the measuring process. Make sure that there are no obstructions, e.g. shelves, above the device.
- Deviations of the room temperature or fluctuations (drafts) lead to changes in the sedimentation values. The device must therefore not be placed in the vicinity of radiators or other heat sources.
- The samples must not be taken directly from the refrigerator and inserted into the Sediplus® S 2000 NX for analysis. To ensure accurate results, it is crucial that the samples are always at room temperature (18 °C to 25 °C or according to country or laboratory specifications).
- To ensure sufficient ventilation for the device's electronics, it must be at least 10 cm away from the wall.
- Keep the measurement optics completely clean at all times. In particular, no blood or plasma splashes may enter the optics.
- Protect the device after use against dust or other contaminants with the dust cover.

- It is essential that direct sunlight or light irradiation be avoided.

Connect the device to the power supply. The mains plug with the power switch and device fuses are located on the left-rear side of the case. Before connecting the device to the mains power, take note of the safety information.

Make sure it is a sufficient distance away from other devices to provide enough operating space and avoid impairing the device's functioning.

You must be able to reach the power switch on the back.

Do not turn on the power supply for the Sediplus® S 2000 NX until you have properly connected all peripheral devices (barcode scanner, PC-AT keyboard or printer) to the device.

The device must be checked for the requirements regarding interference emission and interference immunity in DIN EN 61326-2-6:2013-09 for electromagnetic compatibility (EMC) for *in-vitro* diagnostic devices (IVD).

The electromagnetic environment should be assessed by a specialist before operating this device.

Do not use this device near sources of strong electromagnetic radiation (e.g., unshielded, intentionally operated radio frequency sources) as they can interrupt proper operation.

When connecting and removing the connection of one or more peripheral devices (barcode scanner, PC-AT keyboard, EDP and / or printer), the Sediplus® S 2000 NX must be disconnected from the power supply!

If you notice unexpected operations or damage to programs and/or data, viruses may be present on your computer.

Never use programs or storage media that may contain viruses.
Portable storage media should be checked for viruses with a virus protection program before use.

Use only the S-Sedivette® from SARSTEDT as a patient blood carrier for Sediplus® S 2000 NX.
With other tubes, such as Westergren tubes, you will get incorrect blood sedimentation values.

Be sure to follow the instructions for handling the SARSTEDT S-Sedivette®. Online Instructions for Use: S-Sedivette®; GB 512 at: www.sarstedt.com/en/download/instructions-for-use

Please note that the assessment of measurement results of underfilled S-Sedivettes is at the discretion of the operating personnel, the laboratory management or the attending doctor.

Please also note the interfering factors described in the specialist literature (e.g. paraproteinemia, lipemia, haemolysis), which may have an influence on the measurement of blood sedimentation and, if necessary, check the plausibility of the measurement result in the event of unusual blood samples.

- Only use the S-Sedivettes for the SARSTEDT blood sedimentation system Sediplus® S 2000 NX.
- The surface of the S-Sedivettes must be clean. There must be no residues from labels or powder traces from gloves on it. Impurities produce additional signals which are recorded and can considerably impair determination of sedimentation. The original label or self-adhesive labels at the same height do not interfere with the measurement.
- In addition, it must be ensured that there are no blood drops in the cap during insertion and that the cone in the S-Sedivette® is not wet. This can occur after the blood sample has been taken or after mixing. After starting, the measurement could be influenced by blood dripping down.

- If the determined value of a measurement is higher than that detected during the zero measurement, a negative value is outputted. The measured value will be marked with a minus sign.
- A negative value can only occur due to an incorrect measurement.
- Check the S-Sedivette® for filling level and external damage (dirt, scratches, etc.). Repeat the measurement. If the error persists, please contact the Service Department of DESAGA GmbH in Wiesloch or your local SARSTEDT branch.

If you insert a new S-Sedivette® into the position, the measurement log of the previously removed S-Sedivette® from this position will be immediately deleted.

The data will no longer be available!

It is recommended to use the options for data transmission to an EDP system or a printer, see chapter 2.3.5 or 2.3.6.

The measuring plate must not be obstructed while moving up and down.

There must be no objects on or next to the measuring plate.

If the measuring plate is blocked for any reason, movement is automatically stopped and is only continued when the error message appearing on the display is confirmed.

Be sure to observe the instructions for cleaning the device. Careless cleaning or failure to observe the regulations can lead to malfunctions!

- Disconnect the appliance from the power supply to clean it
- **Exception:** In the menu you will find the menu item Cleaning, see chapter 8.1:
You have the option to move the measuring plate to the highest position to clean the housing parts under the measuring plate and to deactivate the touch display for 30 seconds to clean it.
Do not spray detergent under the measuring plate.
- Clean the housing parts under the measuring plate only with a damp cloth.
- The device is operated with a voltage of 110-230 V.
Do not allow any liquid to enter the instrument during cleaning.
- In principle, only clean your device from the outside!
- Never use abrasive cleaners, aggressive cleaners or solvents.
- Do not allow disinfectants to enter the interior of the device. Do not use disinfectants in spray cans.
- You can remove traces of dust or lint with compressed air.
- Do not open the instrument yourself under any circumstances to remove contamination inside the housing.
- Have the device serviced at regular intervals by the DESAGA Service Department.
- The device must not be autoclaved.

When handling blood as a potentially infectious material, you must always comply with the occupational safety requirements and clean and disinfect the device properly.

Inform yourself of the disinfection measures required for your device in your hygiene department.

In principle, our "*Terms of Delivery and Payment*" apply. These are noted on the back of the invoice.

Warranty and liability claims are excluded if they are attributable to one or more of the following causes:

- Improper use of the device.
- Improper installation, commissioning, operation and maintenance of the device.
- Operation of the device in the event of defective safety devices or improperly installed or inoperable safety and protective devices.
- Failure to observe the instructions for use regarding transport, storage, assembly, commissioning, operation, maintenance, set-up and disposal.
- Unauthorised changes to the device.
- Disasters caused by foreign bodies and force majeure.
- Improperly performed repairs.
- Technical modifications reserved

All serious incidents relating to the product must be reported to the manufacturer and the competent national authority.

1.2 Intended use

Determination of the erythrocyte sedimentation rate* of venous blood samples in the S-Sedivette®. The measurement results are obtained according to the Westergren method.

Usage environment:

The product is intended for use in a professional environment by qualified medical and laboratory personnel.

Functional principle:

The height of the red blood column is determined by scanning an S-Sedivette® (Art. No. 06.1690.xxx) using infrared diodes. The measurement of the erythrocyte sedimentation rate is carried out in a measuring range of 0-65 mm, which corresponds to 0-116 mm according to Westergren. The phase difference ("red blood clot", plasma, air) is clearly recognised by the diodes or the software of the device.

Sample material:

Whole blood anticoagulated with citrate (0.105 mol/l pH 5.5 trisodium citrate /citric acid buffer solution in a mixing ratio of 1:4 (1 part citrate + 4 parts blood) in S-Sedivettes (Art. No. 06.1690.xxx).

*ESR = Erythrocyte sedimentation rate

1.3 EC Declaration of Conformity

The EC Declaration of Conformity can be viewed in various languages on the website of SARSTEDT AG & Co KG www.sarstedt.com; <https://www.sarstedt.com/en/download/declarations-of-conformity/>

2. Introduction

2.1 Description

Based on SARSTEDT's many years of experience in the field of blood sedimentation, an automatic 40-channel blood sedimentation instrument is available through the use of an IR transmission measurement method with the Sediplus® S 2000 NX.

The Sediplus® S 2000 NX has been specially developed for laboratories with centralised blood sedimentation. After the measurements, the results can be fed into to an EDP system. Of course, data can be retrieved manually via the display at any time.

In the optimal combination with the S-Sedivette® blood sedimentation system, blood sedimentation can be carried out automatically and immediately following safe blood collection.

After the mixing process, the measurement is started by inserting the filled S-Sedivette® into a free measuring position.

After you have configured the S-Sedivette® to any measuring station, the device will detect the new S-Sedivette® and the measurement will commence. Likewise, premature removal is detected and displayed as an error.

The display shows the measurement times and results of all positions.

The measuring plate is moved up and down for the measuring process. Each S-Sedivette® is transilluminated by a measuring beam. Behind the S-Sedivette®, the light beam impinges on a detector. The surface of the erythrocyte layer is detected by a change in the light intensity.

Perfect identification of the sample is ensured by the possibility of connecting a barcode scanner. The status display of each of the 40 channels of the device can be called up and is shown on the display.

After completion of the measurement, the result log with 1h and 2h value (optionally also ½ h and 1h value), date, time, ID number and, if necessary, error messages can be output via the display in the display, via the built-in interface to the EDP or via a printer.

2.2 The "blood sedimentation" measurement

'Blood sedimentation, also referred to as erythrocyte sedimentation reaction: ESR = erythrocyte sedimentation rate; can be defined according to CLSI H02-A5 as: distance (length) between the liquid surface of the plasma and the sedimentation surface of the red blood cells (erythrocytes) in the uncoagulated blood in mm.

Physical-chemical principle

Due to the higher density of erythrocytes compared to plasma, they slowly form sediment in anticoagulated blood. However, the sedimentation of erythrocytes is influenced by various factors that can accelerate or slow down the rate of sedimentation.

In addition to the difference in density, the most important influencing variables include erythrocyte aggregation and its electrical charge. The negative charge of the erythrocytes causes them to repel one another and the sedimentation takes place only slowly. A changed composition of the plasma proteins, especially the increase in the concentration of acute phase proteins in inflammatory processes, reduces the zeta potential (repulsive effect of the negative charging) and thereby promotes the aggregation of the erythrocytes, which leads to faster sedimentation behavior of the blood cells and consequently to increased measured ESR values (5). In contrast, the increase in concentration of negatively charged plasma proteins leads to slower sedimentation of the erythrocytes and thus to lower ESR measured values.

The sedimentation should take place at room temperature within a countersinking tube of defined height and defined diameter. Under the influence of gravity, the sedimentation rate is determined after 1 hour and stated in mm-Westergren.

To perform the blood sedimentation, 2.8 ml of venous blood are taken up in the SARSTEDT S-Sedivette® with 0.7 ml of trisodium citrate (0.105 mol/l pH 5.5 trisodium citrate /citric acid buffer solution) (in a mixing ratio of 1:4 (1 part of citrate +4 parts of blood)), thoroughly mixed in the S-Sedivette® and placed vertically in the holder of the Sediplus® S 2000 NX.

The Sediplus® S 2000 NX measures the difference in height between the lower meniscus of the liquid surface of the plasma and the sedimentation surface of the red blood cells. The sedimentation values of the citrate blood are measured at defined points in time, converted into Westergren mm and output after 1 hour in x mm/h according to the standardised Westergren method. In addition, the Sediplus® S 2000 NX also offers the option of outputting the reduction values in x mm/h after ½ and 2 hours.

Reference values should be set at the local level in accordance with the recommendations for reference values. In view of the progressive increase in ESR with age, separate values should be set for each decade of adult life for men and women. Several other variables influence the ESR and can therefore affect the reference values, e.g. haemoglobin levels, medications, menstrual cycle, pregnancy and smoking. The table below lists reference values for the ESRs which can be used as a guide for the establishment of local reference values.

Reference values for the Westergren-ESR method.(1)

ESR mm Westergren / 1h				
Age (years)	Male	Female	upper limit of normal	
			Male	Female
18-30	3.1	5.1	< 7.1	< 10.7
31-40	3.4	5.6	< 7.8	< 11.0
41-50	4.6	6.2	< 10.6	< 13.2
51-60	5.6	9.4	< 12.2	< 18.6
60-70	5.6	9.4	< 12.7	< 20.2
>70	5.6	10.1	< 30	< 35

A normal ESR does not exclude non-inflammatory organ diseases, organ dysfunction and malignant tumours. (4)

- Measurement results must always be checked for plausibility by the operating personnel, the laboratory management or the attending doctor.
- The general regulatory guidelines and statutory regulations for properly determining blood sedimentation must be noted and complied with.

2.3 Setting up the device / Installation

As there are many error sources that can affect the ESR, it is imperative that the information below be noted when setting up the device:

- Vibrations of the filled sedimentation system lead to incorrectly increased sedimentation values. The device must therefore be placed on a fixed table or on a vibration-free surface. For example, the device must not stand on a laboratory table together with a centrifuge.
- The S-Sedivettes may not be scratched, dusty or contaminated under any circumstances.
- Setting up the filled lowering system on an incline leads to increased sedimentation values. The device must therefore be set up in such a way that the S-Sedivettes stand upright and the S-Sedivettes can be loaded from above.
- The measuring plate will move upwards during the measuring process. Maintain sufficient distance to ensure you are able to work unimpeded.
- Deviations of the room temperature or fluctuations (drafts) lead to changes in the sedimentation values. The device must therefore not be placed in the vicinity of radiators or other heat sources.
- The samples must not be taken directly from the refrigerator and inserted into the Sediplus® S 2000 NX for analysis. To ensure accurate results, it is crucial that the samples are always at room temperature (18 °C to 25 °C or according to country or laboratory specifications).
- To ensure sufficient ventilation for the device's electronics, it must be at least 10 cm away from the wall.
- Keep the measurement optics completely clean at all times. In particular, no blood or plasma splashes may enter the optics.
- Protect the device after use against dust or other contaminants with the dust cover.

- It is essential that direct sunlight or light irradiation be avoided.

- Connect the device to the power supply. The mains plug with the power switch and device fuses are located on the left-rear side of the case. Before connecting the device to the mains power, take note of the safety information.
- Make sure it is a sufficient distance away from other devices to provide enough operating space and avoid impairing the device's functioning.
- You must be able to reach the power switch on the back.
- Do not turn on the power supply for the Sediplus® S 2000 NX until you have properly connected all peripheral devices (barcode scanner, PC-AT keyboard or printer) to the device.

- The device must be checked for the requirements regarding interference emission and interference immunity in DIN EN 61326-2-6:2013-09 for electromagnetic compatibility (EMC) for *in-vitro* diagnostic devices (IVD).
- The electromagnetic environment should be assessed by a specialist before operating this device.
- Do not use this device near sources of strong electromagnetic radiation (e.g., unshielded, intentionally operated radio frequency sources) as they can interrupt proper operation.



Figure 1: Rear wall of the Sediplus® S 2000 NX

2.3.1 Technical data:

Manufacturer and device data	
Device:	Sediplus® S 2000 NX
Order no.:	SARSTEDT No. 90.1092
Manufacturer:	SARSTEDT AG & Co. KG
Address:	Sarstedtstraße 1 D-51582 Nümbrecht
Device data	
Measuring principle:	IR transmission measurement method with extraneous light compensation Blood sedimentation measurement in the SARSTEDT S-Sedivette® Conversion to Westergren values
Measurement accuracy:	+/- 1 mm of the measurement path
Measuring range:	0-65 mm, this corresponds to a sedimentation value range of 0-116 mm according to Westergren. Since the S-Sedivette® has a larger diameter and is shorter than the Westergren lowering tube, the determined lowering value is converted into Westergren values and displayed by the Sediplus® S 2000 NX.
Display:	Measuring time in hours and minutes Measured value in millimetres (Westergren values) Touch display for status, information, error messages, as well as the operation of the device 40 x LED status check of the measuring stations
Interfaces:	2x USB for connecting a barcode scanner, PC-AT keyboard or a printer 1 x network socket for transferring measurement logs to EDP
Sample carrier:	SARSTEDT S-Sedivette®
Measuring stations:	40 positions

Data storage:	A built-in battery ensures data retention (measured data, date, time) in the event of a power failure
Data output to:	- Display - Printer - EDP / HL7
Power supply:	110-230 V~, (-10/+10 %), 50-60 Hz
Power consumption:	16 W
Fuses:	0.8 AT (2 fuses, 5 x 20 mm)
Interference emission and immunity to interference	The device must be checked for the requirements regarding interference emission and interference immunity in DIN EN 61326-2-6:2013-09 for electromagnetic compatibility (EMC) for <i>in-vitro</i> diagnostic devices (IVD).
Measurements:	Width: 270 mm Depth: 355 mm Height without S-Sedivettes: 208 mm (retracted measuring plate) height with inserted S-Sedivettes: 305 mm
Weight:	6.7 kg
Permissible environmental conditions when storing the device:	+10 °C to +40 °C, relative humidity max. 80 %, non-condensing 
During operation:	+10°C During operation: Room temperature (18 °C to 25 °C or according to country or laboratory specifications), at relative humidity of max. 80 %, non-condensing.  It is essential that direct sunlight or light irradiation be avoided.

2.3.2 Unpacking, transport and storage

The Sediplus® S 2000 NX (SARSTEDT No. 90.1092) is packed with all accessories in a box.


- Check the integrity of the packaging box and the device.
- In case of transport damage, contact the carrier immediately! Observe the reporting deadline of the individual transport companies (rail, mail, parcel services or freight forwarder). It may be only 24 hours.
- Defects and damages must be reported immediately to SARSTEDT AG & Co. KG!
- Check the scope of delivery based on the list in the next chapter.
- The transport packaging should not be disposed of.

- To transport the device, e.g. to ship it for manufacturer maintenance (see Chapter 7.4 Maintenance at the manufacturer), use only the supplied packaging.
- If the appliance is to be stored in your facility, please observe the environmental conditions specified below.

- It is essential that direct sunlight or light irradiation be avoided.

- Permissible ambient conditions when storing the device: +10° to +40°C, max. 80% rel. humidity non-condensing.

2.3.3 Delivery contents

The Sediplus® S 2000 NX (SARSTEDT No. 90.1092) consists of:

Number	Name
1	Sediplus® S 2000 NX
1	Mains connection cable
1	Dust cover
1	Instructions for Use
1	Test Sedivettes (set of 3)
1	Quick Start Guide

2.3.4 Connecting a PC-AT keyboard and / or a barcode scanner

When connecting and removing the connection of one or more peripheral devices (barcode scanner, PC-AT keyboard, EDP and / or printer), the Sediplus® S 2000 NX must be disconnected from the power supply!

On the back of the Sediplus® S 2000 NX you will find not only the RJ-45 port for creating an Ethernet connection to a network, but also two USB ports for connecting peripheral devices:

You can connect a barcode scanner, a PC-AT keyboard or a printer to the two USB ports.

The PC-AT keyboard and the barcode scanner can be used immediately after switching on the Sediplus® S 2000 NX.

If your barcode scanner does not recognise the submitted barcode, reads incorrectly or attaches characters (for example leading zeros), this is very often due to a different coding system for the barcode being used (for example, transmission of checksums).

In this case, it is necessary to reprogram the barcode scanner to the coding system which you use. Observe the instructions in the enclosed instructions for use of the barcode scanner.

The function of a barcode scanner can be tested in this service menu (menu service barcode scanner). To do this, it must be connected to the USB port, see Figure 2.


Figure 2: Test barcode scanner

2.3.5 Connecting an optional printer

You can connect a barcode scanner, a PC-AT keyboard or a printer to the two USB ports.

Insert the USB connection of the printer into the USB port on the back of the device, see Figure 1. Under the menu item Print settings (see chapter 5.1.7), you have the option of testing the functioning of the printer. You will receive the following printout:

```
Date: 24.06.2022
Time: 09:05:52
Position: 123456
Mode: 0.5h / 1h
Underfilled
ID: TEST PRINT
Measurement 1: 123
Deficient
Delayed
Measurement 2: 456
Deficient
Delayed
```

Figure 3: Test print

2.3.6 Description of the HL7 interface

The Sediplus® S 2000 NX can be connected to a network via the network socket (RJ45) included on the back of the device. If the device is connected to a network via a LAN cable, the IP address of the device is displayed in the service system settings menu, see Chapter 5.2.

To establish a network connection, an integrated DHCP server (Dynamic Host Configuration Protocol) is required, which automatically assigns IP addresses.

The device has a modern interface for transmitting measurement data to a laboratory information system (LIS), which was programmed on the basis of the HL7 standard. This message transmitted by the device is structured as follows:

MSH Segment

Position	Data	Description
0	MSH	Segment name
1		Field separator
2	^~\&	additional separators
3	S2000	Sending application
4	Sarstedt	Sending process
5	{config1}	Receiving Application
6	{config2}	Receiving Facility
7	YYYYMMDDHHMMSS	Time of message creation
8		Blank
9	OML^O21	Message type and event
10	{id}	Message control number
11	P	Processing mode
12	2.6	HL7 Version Number
13		Blank
14		Blank
15	NE	Condition for acknowledgement of receipt
16	NE	Condition for processing confirmation

Example:

```
MSH|^~\&|S2000|Sarstedt|{config1}|{config2}|20201030121501||OML^O21|{id}|P|2.6||NE|NE<CR>
```

PID Segment

Position	Data	Description
0	PID	Segment name
1	1	PID segment number
2		Blank
3	{data1}	Patient ID list, here: Barcode

Example:

```
PID|1||{data1}<CR>
```

ORC segment

Position	Data	Description
0	ORC	Segment name
1	SC	Order control
2	{data1}	Order number, barcode here
3		Blank
4	CM	Order Status

Example:

```
ORC|SC|{data1}|CM<CR>
```

OBR segment

Position	Data	Description
0	OBR	Segment name
1	1 for OBR 1 2 for OBR 2	OBR segment number
2	{data1}	Order number, here: Barcode
3		Blank
4	{config3/4}	Performance identification (test abbreviation)
5		Blank
6		Blank
7	YYYYMMDDHHMMSS	Time of measurement
8		Blank
9		Blank
10		Blank
11		Blank
12		Blank
13		Blank
14		Blank
15		Blank
16		Blank
17		Blank
18		Blank
19		Blank
20	{data2}	Measured value
21	{data3}	If applicable, additional information on the measured value (?)

Example:

```
OBR|1|{data1}||{config3}|||20201030121502|||{data2}|{data3}<CR>
```

A complete HL7 message for the barcode 1234567 with the two measured values 1.11 at 14:35:00 and 2.22 at 15:05:00, measured on 30/10/2020 would then look as follows:

```
<VT>
MSH|^~\&|S2000|Sarstedt|Laboratory 1|Nuembrecht|20201030181500||OML^O21|0001|P|2.6|||NE|NE<CR>
PID|1||1234567<CR>
ORC|SC|1234567|CM<CR>
OBR|1|1234567||Blood sedimentation 30 min|||20201030143500|||1.11|ok<CR>
OBR|2|1234567||Blood sedimentation 60 min|||20201030150500|||2.22|ok<CR>
<FS><CR>
```

"Laboratory 1", "Nuembrecht", "Blood sedimentation 30 min" and "Blood sedimentation 60 min" were assumed for the fields {config 1-4}. A measured value status "ok" was assumed for the field {data3}.

All necessary information to allow integration of the device into your LIS can be entered in the settings menu, see 5.1.8 Network settings / HL7

3. The functional elements of the Sediplus® S 2000 NX


Figure 4: Front view of the Sediplus® S 2000 NX

3.1 The measuring plate

The measuring plate has 40 guide holes for the S-Sedivettes, each with a green control light-emitting diode. Each of the 40 guide holes, every one with a control LED, corresponds to a position for carrying out a measurement.

The S-Sedivettes are inserted vertically through the guide holes into the base of the base module.

The sensor system for sedimentation measurement is accommodated in the measuring plate. Therefore, the measuring plate must not be blocked from moving (moving up and down during the measurement).

Use only the S-Sedivette® from SARSTEDT as a patient blood carrier for Sediplus® S 2000 NX. With other tubes, such as Westergren tubes, you will get incorrect blood sedimentation values.

3.2 The light emitting diodes of the measuring plate

The light-emitting diodes indicate the various states of the associated position or of the entire device.

All 40 light-emitting diodes flash simultaneously.

Optical signalling: The measuring plate is moved up and down in the next 5 seconds for measurement. During this time, no S-Sedivettes may be used. Wait for the measuring process, then insert the S-Sedivettes.

The light-emitting diode of a certain position does not light up and the associated position is unoccupied.

This position is vacant. You can use an S-Sedivette® for measuring.

The light-emitting diode of a position lights up. The position is occupied by an S-Sedivette®.

The measurement with the inserted S-Sedivette® has not yet been completed.

The current status of this position can be queried via the display.

The light-emitting diode of an occupied position does not light up.

The measurement is finished. The value can be retrieved on the display. You can remove use the S-Sedivette® from the device. The measuring position is thus free for the next measurement.

The light-emitting diode of a position lights up, but the position is not occupied.

The S-Sedivette® was removed from this position during the current measurement.

The following warning will appear on the display.


Figure 5: S-Sedivette® removed

- If you have corrected the error using one of the two options, the position is in one of the previously described states.
- If several LEDs have flashed before, you will see the error message of the next position on the display.

3.3 The Display

The Sediplus® S 2000 NX has a modern touch display for operating and displaying the measured values. With the help of the operator panel, you can query all measurement data from completed and ongoing measurements via the display. The device status - date, day of the week, time, number of vacant and filled positions and the number of completed measurements - is displayed. The following figure shows the display of the main view. In the list of measurements you can scroll up and down by touching the touch display, in the lower right corner you can sort the measurement results according to your needs.


Figure 6: Main view

The following information is displayed in the main view:

Table 1

Position	Displays the position in which an S-Sedivette® is located. The position numbering is located on the measuring table cover
ID	Sample ID that uniquely identifies a sample. The ID can be entered in different ways, see, Chapter 4.5 Entering the ID number.
Measurement type	Fixed measurement times
1. value	First measured sedimentation value in mm-Westergren
2. value	Second measured sedimentation value in mm-Westergren
Duration	Duration of measurement. The duration begins as soon as an S-Sedivette® is inserted into the device and only ends when both sedimentation values (1st Value and 2nd value) and/or the S-Sedivette® was removed from the device.

4. Operation of the Sediplus® S 2000 NX

4.1 Switching the device on

Switch on the Sediplus® S 2000 NX at the power switch on the back of the device. You will see the following display view. In the lower right corner of the view, you can see the software version currently installed on the device. You can view these again later under the menu item Service (Chapter 5.2).

Figure 07: Start screen

After a few seconds, the device will switch on automatically. You will see the main view on the display panel.


Figure 08: Home screen with the status display

The device is now ready for operation and you can start a measurement in three different ways, for more details see Chapter 4.3:

1. Scanning a sample ID
2. Inserting an S-Sedivette® into a free slot
3. Pressing the red button **+**

The device has an integrated buffer battery, which guarantees data retention even in the event of unintentional switching off of the device, e.g. in the event of a power failure, see Chapter 6.5.

By clicking on the field in the upper left corner of the display view, you will be taken to the menu view, where you can configure settings and find service functions.


Figure 09: Accessing the menu

The individual menu items are explained in Chapter 5.

4.2 Preparation of samples

Please also note the interfering factors described in the specialist literature (e.g. paraproteinemia, lipemia, haemolysis), which may have an influence on the measurement of blood sedimentation and, if necessary, check the plausibility of the measurement result in the event of unusual blood samples.

- Only use the S-Sedivettes for the SARSTEDT blood sedimentation system Sediplus® S 2000 NX.
- The surface of the S-Sedivettes must be clean. There must be no residues from labels or powder traces from gloves on it. Impurities produce additional signals which are recorded and can considerably impair determination of sedimentation. The original label or self-adhesive labels at the same height do not interfere with the measurement.
- In addition, it must be ensured that there are no blood drops in the cap during insertion and that the cone in the S-Sedivette® is not wet. This can occur after the blood sample has been taken or after mixing. After starting, the measurement could be influenced by blood dripping down.

S-Sedivette® for blood collection

Blood is taken using the SARSTEDT S-Sedivette® with trisodium citrate solution (0.105 mol/l pH 5.5 trisodium citrate /citric acid buffer solution) in a mixing ratio of 1:4 (1 part citrate +4 parts blood) as an anticoagulant.

Observe the Instructions for Use of the SARSTEDT S-Sedivette®

The S-Sedivette® is a closed blood collection system that also serves as a sedimentation tube. You do not need to transfer blood or put it in a separate pipette.

Important:

In principle, each blood sample must be mixed thoroughly without foaming immediately after the blood has been drawn and immediately before it is inserted into the device.

This can be done by hand, so that the mixing ball contained in the S-Sedivette® migrates 5 times through the entire blood sample, or gently and comfortably with a device specially developed for this task by SARSTEDT AG & Co. KG, the Sarmix® M 2000 (see Chapter 9 Ordering Information).

Use the SARSTEDT Sarmix® to mix carefully and gently.

Retention of the blood sample:

Use the obtained blood sample as soon as possible. If the sample is stored at room temperature, the analysis should be carried out within 4 hours. The blood sample must not be clotted.

4.3 Inserting the S-Sedivette® and starting a measurement

You can enter the ID number in different ways (or different sequences), insert the S-Sedivette® and start a measurement.

In principle, the following procedure applies to all types of placement:

- If you have used the S-Sedivettes, the measurement will start automatically.
- 20 seconds after you have inserted the last S-Sedivette®, all LEDs start to flash. No further S-Sedivettes may be used now.
- After a further 5 seconds, the measuring plate scans the newly inserted S-Sedivette®(s).
- In this case, the scanning of the S-Sedivettes begins 25 seconds after the insertion of the last S-Sedivette®.
- **The S-Sedivettes used are in any case scanned at the latest 2 minutes after the insertion of the first S-Sedivette®.**
- After scanning, the display shows the status display with the information on overloaded and free positions and how many measurements have been completed. You can now use new S-Sedivettes and start further measurements.

Option 1:

You work with an optional barcode scanner to enter the patient ID numbers.

- Hold the barcode scanner by the barcode for the S-Sedivette® (for example on the S-Sedivette® or on the patient sheet) without first pressing a button on the touch display and without inserting the S-Sedivette® into an unoccupied position.
- If there is no error at this time (for example, because an S-Sedivette® was removed prematurely), the barcode is automatically detected and the following message appears on the display.


Figure 10: Starting a measurement by scanning an ID

- Put the S-Sedivette® in any vacant position, the measurement will be started automatically
- In this case, you will see the home screen with the status display (Figure 12).
- You can cancel the process by confirming the [Cancel] button as long as the measurement has not been started, i.e. the S-Sedivette® has not been inserted into the device.

Option 2:

The barcode or the patient ID number is available separately (i.e. not only on the S-Sedivette®).

- If the patient ID number is not on the S-Sedivette®, you can insert an S-Sedivette® in a vacant measuring position without first pressing a button.
- Once inserted, the measurement starts automatically. The following display appears:


Figure 11: Starting a measurement by inserting an S-Sedivette®

- Enter the ID number.
You can enter the ID number via the touch display or via an optional PC-AT keyboard and confirm with the [OK] key. Another possibility is to accept the ID number proposed in the display (not consistently assigned to the patient) by pressing the [OK] key.
- The measurement starts automatically and the start screen appears (Figure 12)
- You can cancel the process as long as the measurement has not been started by pressing the [Cancel] key or by removing the S-Sedivette®.

Option 3:

The barcode or patient ID number is on the S-Sedivette®. It must therefore be entered before inserting the S-Sedivette®.

- Press the red button  with the plus sign

Figure 12: Starting a measurement by pressing the "+"key

- The following display appears:


Figure 13: Entering the ID

- Enter the ID number.
You can enter the ID number via the touch display or via an optional PC-AT keyboard and confirm with the [OK] key. Another possibility is to accept the ID number proposed in the display (not consistently assigned to the patient) by pressing the [OK] key.

- The following display appears:


Figure 14: Insert S-Sedivette®

- Place an S-Sedivette® in any vacant position.
- The measurement starts automatically as soon as you have inserted an S-Sedivette® into a vacant slot, then the Start screen with the status display appears (Figure 12)
- You can use the [Cancel] button to cancel the process at any time as long as no S-Sedivette® has been inserted into the device.

4.4 Selecting the measurement times

The device automatically determines the blood sedimentation of a sample after ½ h and after 1 h. You have the option to select the measurement times after 1 h and 2 h. You can choose between selecting this when starting each measurement in the input window (see Figure 15) or setting it globally for all measurements in the settings, see 5.1.3 Measurement settings.


Figure 15: Selecting the measurement times before starting a measurement

4.5 Entering the ID number

As explained in the previous chapter (4.3), there are several ways to enter the ID number. On the one hand, this can be done on request via the touch display, or alternatively via a barcode scanner or a PC-AT keyboard (option).

With the Sediplus® S 2000 NX it is also possible to assign the same ID number several times (this presents a risk of confusion!). In this case, you can only draw distinctions via the position, module or time!

- If you press the red button **+** or insert an S-Sedivette®, the device suggests a consecutive ID number.
- This number, if accepted, will be automatically increased by "one" by the instrument from one measurement to another.
- By pressing the [OK] key, you accept this proposed consecutive number as an ID number (not continuously assigned to the patient).
- You can also enter your own ID number via the touch display or a PC-AT keyboard.
- The entered ID can contain both numbers and letters.

4.5.1 Entering the ID number using a barcode scanner (optional)

Another way to enter the ID number is to scan the number with an optional barcode scanner. (Often the patient ID number is affixed as a barcode on the S-Sedivette® or is available on a patient sheet.)

When reading the barcode from the S-Sedivette®, make sure to hold it vertically. Move the barcode of the S-Sedivette® along the barcode scanner. For information on installing the optional barcode scanner, see chapter 2.3.4 or the user manual of the barcode scanner used.

Take the barcode scanner in your hand and guide it over the barcode on the patient form.


Figure 16: Scan ID Number

4.5.2 Entering the ID number via the touch display

You can enter the ID number directly via the on-screen keyboard. To do this, tap in the ID number field on the display. Then the on-screen keyboard appears and you can enter the desired ID. You can hide the on-screen keyboard via the field in the lower right corner.


Figure 17: On-Screen Keyboard

4.6 Viewing already determined values during a measurement

In the main view, all measurement results are displayed at any time. In general, only measured values converted into Westergren values are displayed.

Depending on the selected option (see chapter 5.1.3 Measurement settings), the measurement value determination for one position is completed after one hour (½h/1h option) or after two hours (1h/2h option). Following this, the two sedimentation values of the position are available.

Position	ID	Meas-type	1st value	2nd value	Duration
1	00000038	0.5h/1h	57mm	112mm	1:14
2	00000039	0.5h/1h	57mm	112mm	1:14
3	00000040	0.5h/1h	57mm	112mm	1:14
4	00000043	0.5h/1h	0mm	0mm	1:14
5	00000044	0.5h/1h	59mm	116mm	1:09
6	00000045	0.5h/1h	0mm	29mm	1:08

Ongoing measurements: 0 | Empty positions: 40
Fri, 24.06.2022 09:48:02

Figure 18: Main view

During a measurement, you will see the status display. You can see how many positions have been filled and how many measurements have been completed.

If you tap on a line on the touch display, a window appears with further details about this measurement. This window appears automatically after the pre-set measurement time has elapsed for a period of 15 seconds, indicating that the S-Sedivette® can be removed from the device (see Figure 19). Completed measurements are highlighted in green in the display view. If a possible error was detected by the device during the measurement, the line of the corresponding measurement is highlighted in red; for details on error detection, see Chapter 6 Faults and operating errors.

Position	ID	Meas-type	1st value	2nd value	Duration
1	00				
2	00				
3	00				
4	00				
5	00				
6	00				

Faulty measurement

Position: Pos 4 ID: 00000043 Measurement type: 0.5h/1h

1st value: 0mm 2nd value: 0mm Duration: 1:30

S-Sedivette® underfilled!

PRINT OK

Ongoing measurements: 0 | Empty positions: 39
Fri, 24.06.2022 10:04:23

Figure 19: Results display

**Important!**

If you insert a new S-Sedivette® into the position, the measurement log of the previously removed S-Sedivette® from this position will be immediately deleted.

The data will no longer be available!

It is recommended to use the options for data transmission to a computer system or a printer, see chapter 2.3.5 or 2.3.6.

4.7 Acoustic support

The Sediplus® S 2000 NX has built-in acoustic support for operation.

A distinction is made between:

- A single tone.
- Two short tones in quick succession.

In principle, the **single tone** is used as a warning signal. For example, in the following cases:

- if you are using an S-Sedivette® and the ID number still needs to be entered;
- if you press the [OK] key after entering the ID number;
- if you use the [OK] key to accept the consecutive ID number;
- if a valid barcode has been detected with the optional barcode scanner.

The **double tone** signals the confirmation of the entered ID number. In addition, the double tone appears in the following cases:

- if you remove an S-Sedivette® before the end of the measurement;
- if the measuring table is prevented from moving

4.8 Reprinting / repeating the printing process

If the printer was activated before the start of the measurements, see 5.1.7 Print settings, the measurement log is automatically printed after the preset measurement time has elapsed.

As long as the S-Sedivette® has not been removed from the device, you have the option of selecting the line with the measurement results on the main view page. In the detail view, you can re-print the log as often as you like by tapping on the print field.


Figure 20: Reprint / Repeat printing

4.9 Deleting measurement results

You have the option to manually delete measurement results from both ongoing and completed measurements.

- You can delete individual ongoing measurements by removing the S-Sedivette® and pressing the "Delete" button on the display, see 6.4 Premature removal of an S-Sedivette®.
- If the S-Sedivette® has already been removed, the measurement log can be called up by tapping the corresponding line in the main view and deleted by pressing the "Delete" button.
- The procedure for deleting all or all completed measurements is described in chapter 5.1.4 Deleting measurements.
- While the S-Sedivette® is on the measuring table, its result cannot be deleted.

5. Menu

Tapping the menu button (see Figure 21) takes you to the menu. Here, you can configure different options and call up service and cleaning functions. In addition, error messages can be displayed and deleted.


Figure 21: Accessing the menu

5.1 Settings

Selecting the menu item "Settings" takes you to the following display view. You can use the touch display to scroll up and down the list of settings.


Figure 22: Settings

5.1.1 Language

You have the option of choosing between the languages German and English

Figure 23: Language settings

5.1.2 Start ID

In the Start ID menu item, you can configure the consecutive ID number to any value. To do this, type in the field with the displayed ID number and enter the desired number, e.g. 0000001, via the display keyboard or the PC-AT keyboard which you may have chosen to connect. The ID number is then incremented by one counter for each newly inserted S-Sedivette®.


Figure 24: Set start ID

5.1.3 Measurement settings

5.1.3.1 Measurement type

In this menu item, you can choose to select the measurement type, i.e. the measurement times after ½ h and 1 h or after 1 h and 2 h. The Sediplus® S 2000 NX offers the possibility to determine and output the ½ h and 1 h measured values instead of the 1 h and 2 h measured values. The device's operation does not change as a result.

Optionally, you have the possibility of querying the measurement type at the start of each new measurement.


Figure 25: Measurement settings

5.1.3.2 "Underfill limit 5 ... 10 mm"

In the menu under Options you will find the menu item "Underfill limit in mm".

This option is used to define the limit in millimetres up to which an underfill in the S-Sedivettes by the Sediplus® S 2000 NX should not be recognised as an error (see also chapter 6.2 Underfilled S-Sedivette®). You can set values between 5 and 10 mm by pressing the "+" and "-" buttons.

Example:

You select the value 8 mm for this option:

- If the S-Sedivette® is only filled up to 8 mm (or less) below the widening of the neck of the S-Sedivette®, the measurement or the measured values will not be marked as underfilled in the result log.
- If the S-Sedivette® is filled up to 9 mm (or more) below the widening of the neck of the S-Sedivette®, the measurement or the measured values are marked as underfilled in the result log.
- The assessment of measurement results of underfilled S-Sedivettes is at the discretion of the operating personnel, the laboratory management or the attending doctor. The reason for this is that sedimentation behavior differs in the case of differently filled S-Sedivettes due to the shifted anticoagulant/blood mixing ratio.
- The S-Sedivettes are then correctly filled when the plunger is fully extended, until it reaches the click position (see Instructions for Use of the S-Sedivette® www.sarstedt.com/en/download/instructions-for-use).
- When breaking off the plunger, make sure that nothing is pushed back into the S-Sedivette®.

Changing the limit is at the discretion of the operating personnel, the laboratory management or the attending doctor. The same applies to the assessment of measurement results of underfilled S-Sedivettes.

When evaluating measured values, note that the maximum evaluation distance extends from the meniscus of the blood column to the upper edge of the label. An extension of the underfill limit therefore shortens the usable measuring range for S-Sedivettes which contain a lower initial volume (see Chapter 6.2 Underfilled S-Sedivette®).

5.1.4 Deleting measurement results

Important!

You cannot undo this function.

When you perform this function, all data of incomplete and completed measurements will be deleted.

You have the option to delete all measurements or all completed measurements. After pressing the corresponding button, a warning message appears, after confirming this, the measurement logs will be deleted.

While the S-Sedivette® is on the measuring table, its result cannot be deleted.


Figure 26: Deleting measurements

5.1.5 Date

The Sediplus® S 2000 NX has a built-in electronic clock.
In the Settings menu, you can set the date, time and the displayed format.


Figure 27: Set date


Figure 28: Set date format

5.1.6 Time

The Sediplus® S 2000 NX has a built-in electronic clock.
In the Settings menu, you can set the date, time and the displayed format.
The date, time and day of the week are already preset in the factory.


Figure 29: Set time

5.1.7 Printer settings

If measurement logs are to be printed via a printer connected to the USB interface, the print interface must be activated in the Settings. The functionality of the printer can be checked in this menu item (test print).


Figure 30: Printer settings

If the printer is activated, the measurement log is automatically printed after a measurement is completed.


Figure 31: Printed measurement log

If printing is not possible, the measurement results will not be erased from the device even after the S-Sedivette® has been removed. After the measuring time has elapsed, a warning message appears, see Chapter 6.4 Premature removal of an S-Sedivette®.


Figure 32: Printing not possible warning message

A warning message appears on the display under the measured values


Figure 33: Results display with printing error

5.1.8 Network settings / HL7

In this menu item, the necessary data is entered to connect the device with an LIS (laboratory information system).

The following data is required:

- Server IP
- Server port
- Received Application
- Received Setup / Process

To set up the interface, please contact your system administrator and the DESAGA service department, see Chapter 7 Service and Maintenance.

If the "Establish connection" function is activated, the measurement log is automatically transmitted after completion of a measurement. If transmission is not possible, an information message appears on the display and the measurement data will not be deleted from the device even after removal of the S-Sedivette®, see Chapter 6.4 Premature removal of an S-Sedivette®.

Figure 34: Network settings HL7


Figure 35: Results display with network error

5.2 Service

The service menu contains menu items that are relevant for the maintenance of the device. These menu items are called up by the manufacturer's service technicians.

5.3 Errors

Possible error messages are displayed in this menu item. These can be deleted as shown.


Figure 36: Delete error message

6. Faults and operating errors

Please also note the interfering factors described in the specialist literature (e.g. paraproteinemia, lipemia, haemolysis), which may have an influence on the measurement of blood sedimentation and, if necessary, check the plausibility of the measurement result in the case of unusual blood samples.

6.1 Fault in the opening and closing of the measuring plate

The measuring plate must not be obstructed while moving up and down.
There must be no objects on or next to the measuring plate.

If the measuring plate is blocked for any reason, movement is automatically stopped and is only continued when the error message appearing on the display is confirmed.

Figure 37: Error message measuring table blocked

Remove the interfering objects and confirm the error message in the display.
The measuring plate moves to the starting position. The measurement is repeated after a few seconds.


6.2 Underfilled S-Sedivette®

The S-Sedivettes must be filled correctly for blood sedimentation. For this purpose, care must be taken to ensure that the SARSTEDT S-Sedivettes are handled correctly during the blood collection (see Instructions for Use of the S-Sedivettes). Online instructions for use: S-Sedivette®; GB 512 at: www.sarstedt.com/en/download/instructions-for-use

- The S-Sedivettes are correctly filled when the plunger is fully extended and when breaking off the plunger rod, nothing is pushed back into the S-Sedivette®.
- In case of underfilling, it is at the discretion of the operator whether or not to call the ESR results into question.
- The blood volume in the S-Sedivette® column influences the sedimentation value due to the shifted anticoagulant/blood mixing ratio.

During the first measurement, the device detects whether the S-Sedivette® is underfilled. The height from which an underfill is detected as an error and documented is configured by the operator in the settings options (Chapter 5.1.3.2 "Underfill limit 5 ... 10 mm"). If the S-Sedivette® is underfilled, this is documented and shown on the display.

Sedimentation values depending on the filling level

P 0:	Level at 20 mbar venous pressure
P 1:	5 mm below P 0 Above this point, an S-Sedivette® is not marked as underfilled
P 2:	10 mm below P 0 Below this point, an S-Sedivette® is always marked as underfilled
Range P 1 to P 2	Possible configuration range for users of the Sediplus® S 2000 NX software for marking an underfilled S-Sedivette®.
P 3:	Top edge label End of the usable measurement window

Distances mm absolute		corresponding to mm Westergren
P 0 to P 3	65 mm	116 mm Westergren
P 1 to P 3	60 mm	99 mm Westergren
P 2 to P 3	55 mm	86 mm Westergren

If the blood sedimentation in the S-Sedivette® is lower than 116 mm Westergren, the device may no longer be able to detect a phase difference. In this case, a value of >116 mm is displayed and, if necessary, the sample can be marked as defective.

A few seconds after the insertion of the S-Sedivette®, the first measurement is carried out. As early as this point in time, the device will detect an underfill of the S-Sedivette®.

If an underfill of an S-Sedivette® has been detected, the measurement is marked in red; this indication remains visible for the entire period of the measurement.


Position	ID	Meas-type	1st value	2nd value	Duration
1	00000038	0.5h/1h	57mm	112mm	1:14
2	00000039	0.5h/1h	57mm	112mm	1:14
3	00000040	0.5h/1h	57mm	112mm	1:14
4	00000043	0.5h/1h	0mm	0mm	1:14
5	00000044	0.5h/1h	59mm	116mm	1:09
6	00000045	0.5h/1h	0mm	29mm	1:08

Ongoing measurements: 0 | Empty positions: 40
Fri, 24.06.2022 09:48:02

Figure 38: Underfilled S-Sedivette®

The indication that an S-Sedivette® is underfilled is also included in the measurement log and displayed for 15 seconds immediately following completion of a measurement.

Position	ID	Meas-type	1st value	2nd value	Duration
1	00				
2	00				
3	00				
4	00				
5	00				
6	00				

Ongoing measurements: 0 | Empty positions: 39
Fri, 24.06.2022 10:04:23

Figure 39: Results display with underfilled S-Sedivette®

6.3 Incorrect measurement

Please note:

If the determined value of a measurement is higher than that detected during the zero measurement, a negative value is output. The measured value will be marked with a minus sign.

A negative value can only occur due to an incorrect measurement.

Check the S-Sedivette® for filling level and external damage (dirt, scratches, etc.). Repeat the measurement. If the error persists, please contact the service department of DESAGA GmbH in Wiesloch or your local SARSTEDT branch.

6.4 Premature removal of an S-Sedivette®

The premature removal of an S-Sedivette® is recognised by the device as an error. This is the case if an S-Sedivette® has been removed from its position, but the measurement has not yet been completed (see also chapter 4.3 Inserting the S-Sedivette® and starting a measurement).

If you remove an S-Sedivette® before the end of the measurement, you will receive a warning on the display.


Figure 40: S-Sedivette® removed

- You cannot start a new measurement until the error has been corrected.
- If desired, discontinue the measurement.
To do this, press the Cancel button. All data of the measurement log (ID number and measured values already determined) will be deleted. The position can be filled once again.
- If you reinsert the S-Sedivette® at the same position, the measurement will be continued.
- The device handles and documents the error in different ways, depending on the time and period for which the S-Sedivette® was removed.
 - a. In the period in which the S-Sedivette® was removed, no measurement was pending.**
 - The measurement is continued.
 - Both measured values can be determined at the correct time.
 - b. The time of a measurement fell in the period in which the S-Sedivette® was removed. No more than five minutes have passed between the time when the measurement should have taken place and the time when the S-Sedivette® was reinserted.**
 - If you insert the S-Sedivette®, the missed measurement will be carried out after 25 seconds.
 - The delay in the measurement is less than 5 minutes and will not be considered or documented further.
 - c. The time of a measurement fell in the period in which the S-Sedivette® was removed. More than 5 minutes have passed between the time at which the measurement should have taken place and the time at which the S-Sedivette® was reinserted.**
 - If you insert the S-Sedivette®, the missed measurement will be carried out after 25 seconds.
 - The delay in the measurement is more than 5 minutes and is documented as follows:
 - Display of the time-related faulty measurement in the list of current measurements:
 - If you remove an S-Sedivette® after completion of the measurement, the sedimentation values for this S-Sedivette® will immediately appear on the display.


Figure 41: Result display with Measurement delayed note

Important!

This log can only be seen for 15 seconds. You can view the measured values of the extracted S-Sedivettes at a later point in time. This is only possible for the values of a position as long as no new S-Sedivette® is used at this position.

Note:

The S-Sedivette® can be removed until a measurement is due without the device interrupting the measurement. This means that if the S-Sedivette® is removed, for example, after the zero measurement, and is set back into its position shortly before the 1-hour measurement, the device catches up with the overdue measurement. The same applies to removal after the 1-hour measurement and timely reset before the 2-hour measurement. The device recovers a measurement as soon as the relevant S-Sedivette®(s) has (have) been reset. The fact that measurements were not taken in time (from 5 minutes measuring time exceeded) is noted in the measurement protocol, but not the duration of the delay.

Since the sedimentation continues, a sedimentation value which is too high will be apparent if the measurement is too late. The user must decide whether a delay for the validation of the measurement results can be accepted. The responsible course of action is to reject the sedimentation values with the late measurement indicated by the Sediplus® S 2000 NX and to carry out a new reduction measurement with a new blood sample.

Note that if the sample is taken between measurements, slight vibrations of the S-Sedivette® can lead to incorrect measurement results.

The case is also classified as premature removal of an S-Sedivette® if the transfer of the measurement log to the EDP system (5.1.8 Network settings / HL7) or to a printer (Chapter 5.1.7 Print settings) has been activated, but no data can be transferred to the EDP or to a printer. In this case, you will also receive a warning message. During this time, the duration of the measurement (in the "Duration" display) will continue to tick up. Measurement times can therefore be achieved which assume values greater than 2 hours (however, the sedimentation values have been determined at the correct time).


Figure 42: Network error (HL7 error) / Print error

6.5 Mains failure

If the device is disconnected from the mains supply during a measurement, no error occurs as long as no measured value is recorded during the period of the mains supply failure.

The built-in battery ensures that the data already determined is preserved (stored). This also applies if you switch off the device after completing all measurements and only want to access the values the next day.

However, if a measured value is detected during the period of the mains failure, the power of the built-in battery is not sufficient to carry out a measurement.

After the mains supply is restored, the device recovers the overdue measurement(s).

If the measurement time is exceeded by more than 5 minutes, the measurement is marked as delayed in the display, on a printout or a transmitted measurement report.

Compare the described problems in chapter 6.4 Premature removal of an S-Sedivette®. The measurement cannot be carried out there because the S-Sedivette® has been removed. The device handles both incorrect operations and errors in the same way.

Example for supply voltage failure:

- You have started a measurement via the touch display or a PC. A few minutes later, the appliance is disconnected from the mains or switched off.
- For example, after 2 hours and 10 minutes, you notice the power interruption and switch the device on again.

In this case, all current measurements will be displayed after switching on the device.

7. Service and maintenance

7.1 Cleaning the device

Two functions can be selected to clean the appliance. If the cleaning position is selected, the measuring table moves to the highest position and remains there until the option is deactivated or the menu item is closed. By fixing the measuring table in this position, the device can be cleaned safely underneath the measuring table. The cleaning position can only be selected if no ongoing measurements are present, since the measurement time provided by the fixing of the measuring table could possibly be delayed.

Furthermore, the touch display can be deactivated for 30 seconds. During this time, the display can be wiped, for example, with a cloth without accidentally activating function and control elements on the touch display.


Figure 43: Cleaning

Be sure to observe the instructions for cleaning the device. Careless cleaning or failure to observe the regulations can lead to malfunctions!

- Disconnect the appliance from the power supply to clean it
- **Exception:** In the menu you will find the menu item Cleaning, see Chapter 7.1:
- You have the option to move the measuring table to the highest position to clean the housing parts under the measuring plate and to deactivate the touch display for 30 seconds to clean it.
- Clean the housing parts under the measuring plate only with a damp cloth. Do not spray detergent under the measuring plate.
- The device is operated with a voltage of 110-230 V. Do not allow any liquid to enter the appliance during cleaning.
- In principle, only clean your device from the outside!
- Never use abrasive cleaners, aggressive cleaners or solvents.
- Do not allow disinfectants to enter the interior of the device. Do not use disinfectants in spray cans.
- You can remove traces of dust or lint with compressed air.
- Do not open the appliance yourself under any circumstances to remove contamination inside the housing.
- Have the device serviced at regular intervals by the DESAGA Service Department (see Chapter 7 Service and Maintenance).
- The device must not be autoclaved.

When handling blood as a potentially infectious material, you must always comply with the occupational safety requirements and clean and disinfect the device properly.

Find out in your hygiene department which disinfection measures are required for your device (e.g. you can fumigate the device in an ethylene oxide chamber or clean it with suitable disinfectants).

Example of a disinfectant:

You can use 70% alcohol for disinfection.

7.2 Quality controls with test Sedivettes

Recommendation:

Have your device serviced by the manufacturer once a year.

Check to see whether your device is functioning correctly at the beginning of each measurement series with the SARSTEDT test Sedivettes. In the event of deviations from the defined sedimentation values, you should have your device checked by the DESAGA GmbH Service Headquarters.

The menu item "Test measurement" (menu – Service test measurement) relates to measurements carried out for the purpose of a functionality test. The menu item can only be selected if no measurements are ongoing.

To carry out a test measurement, you need test Sedivettes, which are included with the device on delivery. You can obtain replacements via the Service Department, see Chapter 7 Service and Maintenance.

- Press the Start measurement button
- Insert the first full test Sedivette into the device
- The device automatically scans the test Sedivette immediately and calculates the offset value
- Remove the test Sedivette and place the second test Sedivette with the middle filling in the same position as the first one.
- The device will scan the test Sedivette again and calculates the first measured value, which is displayed both in millimetres and in mm-Westergren.
- Then remove the second test Sedivette and insert the third test Sedivette with the lowest filling into the same position
- The device will scan the test Sedivette again and calculate the second measured value, which in turn is displayed both in millimetres and in mm-Westergren.


Figure 44: Test measurement

To interpret the results, please compare the values displayed in your device with the following table:

Table 2

Offset:	72-68 mm
Measured value 1:	40 mm - 44 mm
Measured value 1:	56 mm Westergren - 61 mm Westergren
Measured value 2:	63 mm - > 65 mm
Measured value 2:	109 mm Westergren - > 116 mm Westergren

If the values are not within the specified limits, please contact a service technician; for contact information see Chapter 7.4 Maintenance at the manufacturer.

- Measurement results must always be checked for plausibility by the operating personnel, the laboratory management or the attending doctor.
- The general regulatory guidelines and statutory regulations for properly determining blood sedimentation must be noted and complied with.

7.3 Measurement optics

Keep the measurement optics completely clean at all times. In particular, no blood or plasma splashes may enter the optics. Protect the device after use against dust or other contaminants with the dust cover.

7.4 Maintenance at the manufacturer's premises

If you have any questions or problems with the device, please contact the service department of DESAGA in Wiesloch or your local SARSTEDT sales agency.

Do not attempt to open or repair the device.

Important: In most service cases, we require the serial number and the software and hardware version of the control panel and the measuring table electronics. This information can be found in the System sub-point in the service menu, see the Service System Menu (5.2 Service).

Please note that electronic components that contain software are installed both in the control panel (touch display) and in the measuring table.


Figure 45: System

Defective or faulty devices are checked and repaired by our service centre as soon as possible as part of our post-repair service system.

To have this carried out, you should send the defective device with a fault description, the completed contamination questionnaire and a copy of the delivery note in the original packaging to the following address or your local SARSTEDT sales agency.

DESAGA GmbH / SARSTEDT GROUP
Service Department
In den Ziegelwiesen 1-7

D-69168 Wiesloch

Phone: +49 (0) 62 22 / 92 88 65
service.des@desaga-gmbh.de

Recommendation:

Have your device serviced by the manufacturer once a year.

Check to see whether your device is functioning correctly at the beginning of each measurement series with the SARSTEDT test Sedivettes. In the event of deviations from the defined sedimentation values, you should have your device checked by the DESAGA GmbH Service Headquarters.

8. Decommissioning / Disposal

The device, as well as the substances used in connection with this device, must be handled and disposed of properly and professionally in accordance with the applicable legislation. When decommissioning the device, please contact the manufacturer SARSTEDT AG & Co. KG, who will take back the device and dispose of it for you in accordance with professional requirements.

The SARSTEDT S-Sedivette® and blood samples that you use in connection with the Sediplus® S 2000 NX device are possible sources of infection. Make sure you take the necessary disinfection measures and ensure proper and professional disposal!

The Instructions for Use are printed on matt paper (80gm/m²), stapled in the waistband and can be fed into your recycling system.

9. General precautions:

Safety information and warnings

1. General precautions: Use gloves and other general personal protective equipment to protect against blood and possible exposure to pathogens transmissible through biological samples.
2. Handle all biological specimens according to the guidelines and procedures in your facility. In case of direct contact with biological specimens, consult a doctor due to the risk of transmission of HIV, HCV, HBV or other infectious diseases. Always follow the safety guidelines and procedures of your facility

10. Limitations

1. The blood samples should be stored at room temperature. The ESR should be determined within the first 4 hours after collecting the blood.
The specimen can be stored in a refrigerator (4 °C) for a longer period (max. 24 hours). The specimen must then be brought to room temperature before being used.
2. The blood sample must be thoroughly homogenised by careful swirling before measuring the erythrocyte sedimentation rate.
3. The measurement requires an ambient temperature of 18–25 °C and must be protected from vibrations, draughts and direct sunlight or in accordance with country and laboratory-specific specifications.
4. The measuring principle complies with CLSI H2-A5.
5. The measuring accuracy is +/-1 mm of the measuring path.
6. The mixing ratio of 1:4 has a direct effect on the analysis result and must be adhered to.
7. Use only with the S-Sedivette® (see note on page 7).

11. Order information

Device / Accessories	Order No.:
Sediplus® S 2000 NX 230V	90.1092
Barcode scanner for Sediplus® S 2000 NX 230V	90.1092.730
Thermal printer for Sediplus®	90.1092.720
Paper roll for thermal printer, 5 pieces	90.188.055
Sarmix® M 2000, 230 V, incl. power supply	90.180.600
Block rotor for 40 S-Sedivettes	92.180.615
S-Sedivette®, EU colour code	06.1690.001
S-Sedivette®, ISO colour code	06.1690.100

Further BSG system solutions are available from SARSTEDT AG & Co. KG.

Article:	Order No.:
S-Monovette® ESR, 2 ml, 66 x 11 mm, with transparent label, EU colour code	05.1079
S-Monovette® ESR, 2 ml, 66 x 11 mm, with paper label, EU colour code	05.1079.001
S-Monovette® ESR, 2 ml, 66 x 11 mm, with paper label, ISO colour code	05.1079.100
Plastic countersunk pipette with O-marking	86.1996
Plastic countersink pipette with printed scale	86.1996.062
ESR rack for the S-Monovette® ESR, with graduated rear wall	90.1060
ESR rack for the S-Monovette® ESR, without rear wall	90.1060.062
Microvette® CB 200 ESR, EU colour code	18.1325
Microvette® CB 200 ESR, ISO colour code	18.1325.100
ESR rack for the Microvette® CB 200 ESR with graduated real wall	90.1091

If you have any questions about other SARSTEDT products or the SARSTEDT product range, please contact:

SARSTEDT AG & Co. KG
Sarstedtstraße 1
D-51582 Nümbrecht

Phone: +49 2293 305 0
Fax: +49 2293 305 3450

info@sarstedt.com
www.sarstedt.com

or your local SARSTEDT sales agent.

12. Warranty and liability

In principle, our "Terms of Delivery and Payment" apply. These are noted on the back of the invoice. Warranty and liability claims are excluded if they are attributable to one or more of the following causes:

- Improper use of the device.
- Improper installation, commissioning, operation and maintenance of the device.
- Operation of the device in the event of defective safety devices or improperly installed or inoperable safety and protective devices.
- Failure to observe the instructions for use regarding transport, storage, assembly, commissioning, operation, maintenance, set-up and disposal.
- Unauthorised changes to the device.
- Disasters caused by foreign bodies and force majeure.
- Improperly performed repairs.
- Technical modifications reserved

All serious incidents relating to the product must be reported to the manufacturer and the competent national authority.

13. Glossary

Anticoagulant	Trisodium citrate solution (0.105 mol/l pH 5.5 trisodium citrate /citric acid buffer solution, equivalent to 3.2% trisodium citrate / citric acid buffer solution in a mixing ratio of 1:4 (1 part citrate + 4 parts blood), 0.7 ml predosed
Barcode	Patient ID, number on S-Sedivette® or patient form
Barcode scanner	Scanner for recording ID numbers
Erythrocyte sedimentation rate (ESR)	Erythrocyte sedimentation rate in x mm/h, in short: blood sedimentation
Documentation	Analysis results with patient name and temperature
Measuring stations	40
Identification number	Device-specific distribution of the analysis number, any number entered or patient barcode
Measuring plate	Guide plate for 40 S-Sedivettes with built-in sensors and control light-emitting diodes
Measurement log	Results with 1/2 h value (½/1 h), date, time and ID number Output via touch display, serial printer or to the EDP
Measuring points	40 Specimen receptacles for S-Sedivette®
Room temperature	Temperature range, 18-25°C according to CLSI H02 – A5(1), or according to country or laboratory-specific specification)
HL7 interface	Interface for transmitting measurement data to a laboratory information system (LIS)
Sedimentation speed	Sedimentation of erythrocytes in x mm/h
S-Sedivette®	SARSTEDT secure blood collection system
Westergren	Sedimentation value in mm/h

14. List of illustrations

Figure 1: Rear wall of the Sediplus® S 2000 NX	67
Figure 2: Test barcode scanner.....	70
Figure 3: Test print.....	71
Figure 4: Front view of the Sediplus® S 2000 NX.....	74
Figure 5: S-Sedivette® removed.....	75
Figure 6: Main view.....	76
Figure 07: Start screen	77
Figure 08: Home screen with the status display	77
Figure 09: Accessing the menu.....	78
Figure 10: Starting a measurement by scanning an ID	79
Figure 11: Starting a measurement by inserting an S-Sedivette®	80
Figure 12: Starting a measurement by pressing the "+"key	80
Figure 13: Entering the ID	81
Figure 14: Insert S-Sedivette®	81
Figure 15: Selecting the measurement times before starting a measurement	82
Figure 17: On-Screen Keyboard	83
Figure 16: Scan ID Number	83
Figure 18: Main view	84
Figure 19: Results display	84
Figure 20: Reprint / Repeat printing	86
Figure 21: Accessing the menu.....	87
Figure 22: Settings	87
Figure 23: Language settings.....	87
Figure 24: Set start ID.....	88
Figure 25: Measurement settings.....	88
Figure 26: Deleting measurements.....	89
Figure 27: Set date	90
Figure 28: Set date format	90
Figure 29: Set time	90
Figure 30: Printer settings	91
Figure 31: Printed measurement log	91
Figure 32: Printing not possible warning message	91
Figure 33: Results display with printing error	92
Figure 34: Network settings HL7	92
Figure 35: Results display with network error.....	92
Figure 36: Delete error message	93
Figure 37: Error message measuring table blocked.....	93
Figure 38: Underfilled S-Sedivette®	95
Figure 39: Results display with underfilled S-Sedivette®	95
Figure 40: S-Sedivette® removed.....	96
Figure 41: Result display with Measurement delayed note	97
Figure 42: Network error (HL7 error) / Print error.....	97
Figure 43: Cleaning.....	98
Figure 44: Test measurement.....	100
Figure 45: System.....	101

15. Literature

- 1) CLSI* H2-A5 2011: Procedures for Erythrocyte Sedimentation Rate Test.
- 2) Thomas, L.: Labor und Diagnose. Digitale Ausgabe Labor & Diagnose 2020; Kapitel 19.3.2 Blutkörperchensenkungs-Reaktion (BSR) (www.labor-und-diagnose-2020.de/k19.html#_idTextAnchor8638)
- 3) ICSH recommendations for measurement of erythrocyte sedimentation rate. J. Clin. Pathol. 1993 46:198-203
**CLSI (Clinical and Laboratory Standards Institute)*
- 4) Bruhn et al.; Lehrbuch der Labormedizin; [Textbook of Laboratory Medicine] Schattauer, 1999; ISBN 3-7945-1856-X
- 5) Gressner et al; Lexikon der Medizinischen Laboratoriumsdiagnostik [Encyclopaedia of Medical Laboratory Diagnostics] Reference work; 3rd Edition. 2019; Springer. ISBN 978-3-662-48987-1

16. Key for symbols and labels:

These Instructions for Use contain the following important symbols to indicate hazards and operating errors. The chapter "Safety information" contains detailed information on safety when handling the Sediplus® S 2000 NX.

This "Important" sign means that in the section:

- Important notes are provided.
- An immediate or potentially imminent danger to the life and health of persons is described.
- An indication of dangerous situations can be found.

Failure to observe this symbol may result in minor to serious or fatal injury or damage to property.

Warning of dangerous electrical voltage!

This "info symbol" indicates information on the proper use of the device or general information to provide a better understanding. Failure to observe the instructions may result in the device being operated incorrectly or even being damaged.

This "tip symbol" provides you with tips on use or particularly useful information on how to use the device in the optimal manner.

Follow the instructions for use

Permissible temperature range

Keep away from sunlight

Store in a dry place

Biological hazard

Article number

Batch number

CE mark

In vitro diagnostic device

Manufacturer

Country of manufacture

Date of manufacture

Unique product identification

Serial number

Separate collection of electrical and electronic equipment

All serious incidents relating to the product must be reported to the manufacturer and the competent national authority.

17. Contamination questionnaire

Before you send the Sediplus® S 2000 NX device to DESAGA GmbH /SARSTEDT group for repair, please complete the contamination questionnaire.

DESAGA GmbH / SARSTEDT-GRUPPE
Service Department
In den Ziegelwiesen 1-7
D-69168 Wiesloch
Tel.: +49 (0) 62 22 / 92 88 65
service.des@desaga-gmbh.de

Send the device together with the completed contamination questionnaire, a brief description of the problem encountered and a copy of the delivery note to the service department of DESAGA GmbH / SARSTEDT GROUP Wiesloch or to your local SARSTEDT sales agency.

Table 3: Contamination questionnaire

Contamination questionnaire for repair orders	
Dear Customer,	
for safety reasons, we ask you to answer all the following questions before submitting your device to DESAGA GmbH / SARSTEDT GROUP Wiesloch for maintenance, repair or return.	
Company:	Department:
City:	Street:
Name:	Phone No.:
Device/Item:	SN:
<input type="checkbox"/> We hereby confirm: The device is not contaminated. The device is free of harmful substances	
<input type="checkbox"/> The device is free of harmful substances	
Substance class	Substance name
<input type="checkbox"/> Toxic substances	
<input type="checkbox"/> Corrosive substances	
<input type="checkbox"/> Explosive substances	
<input type="checkbox"/> Radioactive material	
<input type="checkbox"/> Contagious substances	
<input type="checkbox"/> Flammable substances	
<input type="checkbox"/> Other dangerous substances	
<input type="checkbox"/> The device has been decontaminated in accordance with the legal regulations	
<input checked="" type="checkbox"/> Please tick the appropriate box!	
Description of the decontamination carried out:	
Date:	Signature:

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