


---

## For the use of Sarstedt consumables in cryopreservation

---

Cryo Performance Tested items fulfil the following certified purity criteria:

- ✓ Sterile
- ✓ Pyrogen-free/endotoxin-free
- ✓ Non-cytotoxic
- ✓ Non-mutagenic
- ✓ CE 

We guarantee that the following limit values are met:

Sterility validated according to ISO 11137

Pyrogens/endotoxins <0.06 EU/ml

Non-cytotoxic according to ISO 10993

Non-mutagenic according to the Ames Test II