

Quantidade de poços	96 pré-carregados	96	96	96	96	384	96	384
Meia saia/Saia	sem	sem	Meia	Completo	Meia	Completo	Meia	Completo
Perfil	High	High	High	Low	Low	-	"Lightcycler"	"Lightcycler"
Ref. ^a placas para PCR	72.985	72.1978 72.1978.010	72.1979 72.1979.010 72.1979.003 72.1979.201 72.1979.700 72.1979.102 72.1979.132	72.1980 72.1980.010 72.1980.201 72.1980.600 72.1980.700	72.1981 72.1981.010	72.1984.202	72.1982.202	72.1985.202

BioRad®/MJ Research®

CFX96 Touch™ Real-Time PCR				●				
CFX384 Touch™ Real-Time PCR						●		
CFX Automation System II								
T100™ Thermal Cycler	●	●		●	●			
S1000™ Thermal Cycler		●	●	●		●		
C1000Touch™ Thermal Cycler		●	●	●	●	●		
iCycler iQ™ Thermal Cycler	●	●	●					
iQ4™ Thermal Cycler	●	●	●					
iQ5™ Thermal Cycler	●	●	●					
MyCycler™ Thermal Cycler	●	●	●					
Chromo4™		●		●				
Opticon™, Opticon2™				●				
BaseStation™				●				

Corbett Research®/Qiagen®

Palm Cycler 96 Well		●						
Palm Cycler 384 Well				●				

Eppendorf®

Mastercycler® nexus	●	●						
Mastercycler® ep realplex		●	●			●		
Mastercycler® gradient	●	●		●				
Mastercycler® ep gradient	●	●	●	●				
Mastercycler® pro	●	●		●				

Ericom®

Deltacycler			●	●				
SingleBlock			●	●				
TwinBlock			●	●				

MWG®

Primus 96 Well		●		●				
Primus 384 Well						●		
The Q-Lifecycler		●	●	●				

Roche®

Lightcycler® 96 System							●	
Lightcycler® 480 System							●	●

Stratagene®/Agilent®

AriaMx Real-Time PCR System				●	●			
Mx3000P™	●	●						
Mx3005P™	●	●	●	●				
Mx4000™	●	●	●					
Gradient Cycler		●		●				
Robocycler 384 Well							●	

Techné®

Cyclogene		●		●				
Flexigene		●	●	●				
Genius/Genius Quad		●	●	●				
OMN-E		●	●	●				
PCR Express	●	●	●				●	
Primus 96		●						
Px2/PxE		●	●				●	
Quantica			●	●				
TC412/TC512		●		●			●	
Touchgene/Touchgene Gradient		●	●	●	●			